ottobock.

C1000 DS

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EN Instructions for Use

Zusatzoptionen C1000 DS

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Additional options for the C1000 DS

The power wheelchair can be equipped with the following additional options (applicable options are checked):

| Electronic drive-away lock* | | | |
|---|-------------------------------|--------------|--|
| [] Function enabled | [] Function disabled | | |
| If enabled, the function is activated by pressing the mode key on the control panel. | | | |
| The function is deactivated with the joystick. | | | |
| | | | |
| Electronic steering lock* | | | |
| [] Function enabled | [] Function disabled | | |
| | | | |
| 2-way adapter cable* | | | |
| [] Adapter cable included | [] Adapter cable not included | | |
| The included adapter cable allows the following functions to be selected directly through separate buttons: | | | |
| [] Footrest, left | [] Footrest, right | [] Seat tilt | |
| [] Seat height adjustment | [] Back angle adjustment | [] Light | |
| [] Direction indicator, right | [] Direction indicator, left | [] Horn | |
| [] Warning flasher | [] other: | | |
| | | | |

^{*}See the "Accessories" chapter for more information.

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Instructions for Use

Table of contents

| 1 | General Information | 117 |
|-----|---|-----|
| 1.1 | Foreword | 117 |
| 1.2 | Intended use | 117 |
| 1.3 | Field of application | 118 |
| 1.4 | Service | 118 |
| 2 | Safety instructions | 119 |
| 2.1 | Explanation of Symbols | 119 |
| 2.2 | Standards and directives | 120 |
| 2.3 | General safety information | 120 |
| 2.4 | Safety requirements for transportation, | |
| | storage and assembly | |
| 2.5 | Safety requirements for operation | 124 |
| 2.6 | Safety requirements for care, maintenance | |
| | and disposal | 130 |
| 2.7 | Requirements for the user | 132 |
| 2.8 | Safety functions | 132 |
| 2.9 | Warning symbols and type plates | 134 |
| 3 | Product description | 136 |
| 4 | Delivery and preparation for use | 137 |
| | | |

| 4.1 | Delivery | . 137 |
|-------|---|-------|
| 4.2 | Settings | . 138 |
| 4.3 | Initial operation | . 138 |
| 5 | Transportation and storage | 140 |
| 5.1 | Transportation in a Wheelchair Accessible Vehicle | . 140 |
| 5.1.1 | Required Accessories | . 141 |
| 5.2 | Transport of the Power Wheelchair | . 141 |
| 5.2.1 | Preparing for transportation | . 142 |
| 5.3 | Storage | . 143 |
| 6 | Operation | . 144 |
| 6.1 | Adjustment possibilities | . 144 |
| 6.1.1 | Back rest | . 145 |
| 6.1.2 | Adjusting the luggage carrier | . 145 |
| 6.1.3 | Side panel with arm rest | . 146 |
| 6.1.4 | Control panel | . 148 |
| 6.1.5 | Footrests | . 149 |
| 6.2 | Getting into and out of the power wheelchair | . 150 |
| 6.2.1 | From/to the side | 150 |
| | 1 10111/10 1116 3146 | |

| Control unit | 52 |
|---|---------------|
| Control panel 1 | 52 |
| Separate LCD monitor 1 | 54 |
| Switching on and off | 54 |
| Driving function 1 | 55 |
| "Battery Capacity" indicator | 57 |
| Controlling additional electrical functions 1 | 58 |
| Drive-away lock1 | 59 |
| Steering Lock | 61 |
| Lighting 1 | 62 |
| Releasing and engaging the wheel lock 1 | 63 |
| Batteries 1 | 65 |
| Charging the batteries 1 | 66 |
| Charger 1 | 67 |
| Removing the battery cover 1 | 69 |
| Replacing the batteries1 | 70 |
| Changing the battery installation position 1 | 72 |
| Accessories1 | 74 |
| Seat options1 | 74 |
| Seat height adjustment | 75 |
| | Control panel |

| 7.1.2 | Electric seat tilt | 178 |
|--------|---|-----|
| 7.1.3 | Electric back angle adjustment | 178 |
| 7.1.4 | Recaro® seat | 179 |
| 7.1.5 | Contour seat | 181 |
| 7.2 | Push button module | 181 |
| 7.3 | 2x adapter cable for Piko or Buddy button | 183 |
| 7.4 | Mechanically elevating footrests | 184 |
| 7.5 | Electrically adjustable footrests | 184 |
| 7.6 | Control panel holder | 185 |
| 7.6.1 | Swing-away control panel holder with | |
| | removable control panel | 185 |
| 7.6.2 | Height-adjustable control panel holder | 185 |
| 7.7 | Special controls | 186 |
| 7.8 | Attendant control | 186 |
| 7.8.1 | Functional overview | 187 |
| 7.9 | Additional indicator and control elements | 189 |
| 7.9.1 | Separate LCD monitor with infrared | 189 |
| 7.9.2 | External kilometre indicator | 192 |
| 7.10 | Additional options | 192 |
| 7.10.1 | 1 Lap belt | 192 |
| | | |

| 7.10.2Four-way chest strap/chest belt (for Recaro® seats only) | | | |
|---|---|--|--|
| 7 10 3 | - | | |
| 7.10.0 | 7.10.3Headrest mounting adapter (standard seat, contour seat) | | |
| 7.10.4 | Other optional add-on components | | |
| 8 | Malfunctions / troubleshooting 195 | | |
| 8.1 | Warning 196 | | |
| 8.2 | Error | | |
| 8.2.1 | Error Indication: Control Panel / LCD Monitor 197 | | |
| 8.2.2 | Error Indication: Attendant Control | | |
| 8.3 | Defect/failure | | |
| 9 | Maintenance, cleaning and disinfection 202 | | |
| 9.1 | Maintenance intervals | | |
| 9.2 | Changing the fuse | | |
| 9.3 | Wheel replacement | | |
| 9.3.1 | Steering caster replacement | | |
| 9.3.2 | Drive wheel replacement | | |
| 9.3.3 Replacing the casing and/or tube for | | | |
| | pneumatic tyres | | |
| 9.4 | Changing defective lights | | |
| 9.4.1 | Front lights | | |

| 9.4.2 | Rear lights | 210 |
|-------|-----------------------|-----|
| 9.5 | Cleaning and Care | 211 |
| 10 | Technical data | 212 |
| 11 | Disposal | 215 |
| 12 | Information on re-use | 215 |
| 13 | Liability | 216 |
| 14 | CE Conformity | 216 |
| | | |

List of Figures

| Figure 1 | Main components |
|-----------|--|
| Figure 2 | Anti-tipper |
| Figure 3 | Fuse |
| Figure 4 | Inserting the cross bolt into the holder; locked cross bolt |
| Figure 5 | Eyebolts front, eyebolts rear 142 |
| Figure 6 | Folding the back rest143 |
| Figure 7 | Folding size |
| Figure 8 | Release strap for the back rest 145 |
| Figure 9 | Adjusting the luggage carrier146 |
| Figure 10 | Loosening the thumb screw for removal of the side panel with arm rest146 |
| Figure 11 | Adapting the arm rest to the upper arm length147 |
| Figure 12 | Adapting the arm rests to the forearm length 147 |
| Figure 13 | Adapting the control panel148 |
| Figure 14 | Footrest locking mechanism 149 |
| Figure 15 | Adjusting the lower leg length 150 |
| Figure 16 | Footrests folded up |
| | |

| Figure 17 | Footrests removed | 151 |
|-----------|---|-----|
| Figure 18 | Control panel | 152 |
| Figure 19 | Control panel connections on the underside | 153 |
| Figure 20 | LCD monitor with all symbols | 154 |
| Figure 21 | Information window with drive-away lock | 160 |
| Figure 22 | Front lights | 162 |
| Figure 23 | Rear lights, reflectors | 163 |
| Figure 24 | Releasing the wheel lock | 164 |
| Figure 25 | Charging receptacle | 169 |
| Figure 26 | Removing the battery cover | 170 |
| Figure 27 | Unlocking the battery drawer and removing the track | 171 |
| Figure 28 | Removing the batteries by the handles | |
| Figure 29 | Battery stop | 172 |
| Figure 30 | Release bolt engaged | 172 |
| Figure 31 | Battery stop | 173 |
| Figure 32 | Creep speed indicator on the control unit display / LCD display | 177 |
| Figure 33 | Seat with activated seat height adjustment function | 178 |
| Figure 34 | Electric seat tilt | 178 |

C1000 DS

| Figure 35 | Electric back angle adjustment | 179 |
|-----------|---|-----|
| Figure 36 | C1000 DS power wheelchair with | |
| | Recaro® seat | 179 |
| Figure 37 | Illustration of the push buttons | 180 |
| Figure 38 | Contour seat | 181 |
| Figure 39 | Push button module | 183 |
| Figure 40 | Mechanically elevated footrest | 184 |
| Figure 41 | Elevated footrest | 185 |
| Figure 42 | Control panel holder, swing-away | 185 |
| Figure 43 | Control panel holder, height-adjustable | 186 |
| Fig. 44 | Attendant control | 186 |
| Fig. 45 | Attendant control overview | 187 |
| Figure 46 | LCD module with infrared | 190 |
| Figure 47 | Kilometre indicator | 192 |
| Figure 48 | Adjusting the lap belt | 193 |
| Figure 49 | Adapter for headrest installation kit | 194 |
| Figure 49 | Fuse holder | 205 |
| Figure 50 | Jack | 206 |
| Figure 51 | Removing the steering caster | 207 |
| Figure 52 | Removing the drive wheel | 208 |
| | | |

| Figure 53 | Removing the casing | 209 |
|-----------|--|-----|
| Figure 54 | Attaching the rear light to the steering | |
| | caster fork | 210 |
| Figure 55 | Removing the rear light | 210 |
| Figure 56 | Rear light connector plug | 211 |

1 General Information

INFORMATION

Date of the last update: 2012-08-15

- Please read this document carefully.
- Follow the safety instructions.

INFORMATION

The present Instructions for Use can be viewed at www.ot-tobock.com and downloaded from there. It is possible to increase the display size of the PDF document stored there. For other questions concerning the Instructions for Use please contact the specialist who delivered the product to you.

1.1 Foreword

These instructions for use provide the user as well as his or her attendants with all the required knowledge on the design, function, operation and maintenance of the C1000 DS power wheelchair from Otto Bock Mobility Solutions GmbH. The instructions contain all the information that is needed to use the power wheelchair safely, to determine the possible causes of a malfunction and to help eliminate it.

Knowledge of these instructions for use is absolutely necessary for the safe operation of the power wheelchair. Therefore, the user and his or her attendants must read the instructions for use thoroughly, especially the section on "Safety", before using the power wheelchair. This will ensure that all features of the power wheelchair can be utilised to their full advantage.

1.2 Intended use

The C1000 DS power wheelchair is intended exclusively for individual indoor and outdoor self-transportation by persons with walking impediments or walking inabilities. The C1000 DS power wheelchair may only be combined with the options mentioned in these instructions for use. Otto Bock Mobility Solutions GmbH assumes no liability for combinations with medical devices and/or accessories from other manufacturers outside of the modular system. The C1000 DS power wheelchair may be used for transportation by an attendant if the corresponding optional special controls are installed.

Any other use is considered improper use. The manufacturer is not liable for any personal injury or damage to property resulting from improper use; in such cases, the user has sole liability.

The C1000 DS power wheelchair may only be used by properly trained persons. Training in the use of the C1000 DS power wheelchair is required to protect persons from danger and to ensure that the C1000 DS power wheelchair is operated safely and correctly.

Safe operation of the C1000 DS power wheelchair can only be assured if it is used for the intended purpose in accordance with the information in these instructions for use. The user is ultimately responsible for accident-free operation.

1.3 Field of application

The versatility of the C1000 DS power wheelchair and its modular design make it suitable for users who have walking impediments or walking inabilities due to:

- Paralysis
- Loss of limbs
- Limb defects / deformities
- Joint contractures / defects
- Other disorders

The C1000 DS power wheelchair was designed especially for users who are able to move independently in such a power wheelchair.

The following points should also be considered for individual fittings:

- Body height and body weight (standard version: max. load capacity 140 kg (309 lbs) / optional 200 kg (440 lbs); version with seat height adjustment: 130 kg (287 lbs) / optional 200 kg (440 lbs); version with higher-performance batteries 94 Ah (C5): max. load capacity 140 kg (309 lbs), when using the seat height adjustment 130 kg (287 lbs))
- Physical and psychological limitations
- Age of the patient
- Living conditions
- Environment

1.4 Service

INFORMATION

Service and repairs on the C1000 DS power wheelchair must only be carried out by specialist dealers with personnel authorised and trained by Ottobock. In case of problems, contact the specialist dealer who adapted your C1000 DS power wheelchair for you.

Should you have any questions or a problem that cannot be

resolved despite consulting the instructions for use, please contact Ottobock (see back side of cover for address).

Ottobock endeavours to support their customers in all respects in order to keep them satisfied with their product for a long time.

2 Safety instructions

2.1 Explanation of Symbols

⚠ WARNING

Warnings regarding possible risks of severe accident or injury.

⚠ CAUTION

Warnings regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

INFORMATION

Information regarding operation. Information for the user and service personnel.

2.2 Standards and directives

All safety information contained in these instructions for use refers to the currently valid national laws and regulations of the European Union. In other countries, compliance with the applicable laws and national regulations is required.

In addition to the safety instructions contained in these instructions for use, the user is required to observe and comply with the generally applicable regulations of the BGV (Employer's Liability Insurance Association), UVV (accident prevention regulations) and environmental protection regulations. All information contained in these instructions for use must be complied with at all times without restrictions.

The C1000 DS power wheelchair has been constructed in accordance with state-of-the-art technology and is safe to operate. The safety of the C1000 DS power wheelchair is confirmed by the CE symbol and the declaration of conformity.

2.3 General safety information

△ WARNING

Risk of suffocation. Packaging materials must be kept out of reach of children.

⚠ CAUTION

Risk of accidents and injury if safety instructions are not observed or complied with.

All safety instructions contained in these instructions for use and in all other applicable documents must be observed and complied with. The instructions for use must be available to the user at all times.

⚠ CAUTION

Risk of accidents and injury due to improper use. The power wheelchair must only be used in accordance with the guidelines for proper use. The power wheelchair may only be used by properly trained persons. Only a single person may be transported with the power wheelchair at one time.

⚠ CAUTION

Risk of burns in the proximity of fire. The back upholstery and seat cushion are flame retardant but nevertheless flammable. Therefore utmost CAUTION is required near any sources of open flame or sparks, especially lit cigarettes.

⚠ CAUTION

Risk of hypothermia or burns on wheelchair parts.

Parts can become extremely hot or cold due to extreme temperatures.

Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold) in order to prevent injuries by coming into contact with the parts.

INFORMATION

Only use original options supplied by the manufacturer. The optional components may only be mounted as described here. Failure to comply will void the warranty.

INFORMATION

The serial numbers required for inquiries and orders are found on the nameplate (see Section 2.9).

2.4 Safety requirements for transportation, storage and assembly

Transportation and storage

⚠ WARNING

Risk of accident and injury due to improper use for transportation in wheelchair accessible vehicles. The power wheelchair may only be used for transportation in a wheelchair accessible vehicle if the safety components (such as lap belt) which are offered by Ottobock and suitable restraint systems are used.

Only one person may be transported in the C1000DS at any one time. For optimum protection of the occupants in case of an accident, use the seats and corresponding restraint systems installed in a wheelchair accessible vehicle while it is in operation.

For more information on the use of the C1000DS as a seat in a wheelchair accessible vehicle, please refer to our brochure "Using your Wheelchair / Mobility Base with Seating Shell or Buggy for Transportation in Wheelchair Accessible Vehicles", order number 646D158.

⚠ WARNING

Risk of accident and injury due to incorrect transportation in aeroplanes. In aeroplanes, the power wheelchair has to be transported in compliance with IATA regulations (International Air Transport Association). To this end, the fuse must always be removed and the battery connectors must be insulated so that they are short-circuit-proof prior to check-in of the power wheelchair as luggage.

Note that those batteries in particular which may leak or will not be transported upright must be removed and packaged so they cannot leak or short circuit.

For more information please visit the www.iata.org website. Ottobock recommends contacting the airline directly before every flight to obtain information regarding special transport regulations.

NOTICE

Risk of damage due to improper transportation. Hoisting devices used for transportation must have a sufficient capacity. The power wheelchair must be secured in accordance with the requirements for the transport aid. Only attach tensioning straps to the corresponding eyebolts.

During transportation on lifting platforms or in lifts, the control unit of the power wheelchair must be turned off. Engage the brake. Ensure that the power wheelchair is centred on the lifting platform. None of the power wheelchair components or parts are allowed to be in the danger zone.

NOTICE

Battery damage due to deep discharge. Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

INFORMATION

The tyres of the power wheelchair contain chemical substances that may react with other chemical substances such as cleaning agents or acids.

Black tyres contain soot particles that can cause discolouration and black marks. A suitable surface is required when parking the wheelchair for extended periods of time.

Assembly

⚠ CAUTION

Risk of accidents due to unsecured screw connections. When screw connections with thread locks are loosened, they must be replaced with new ones or secured with a medium-strength liquid thread lock substance (e.g. Loctite 241®).

After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened. Observe the applicable torque specifications.

⚠ CAUTION

Risk of tipping due to an incorrectly mounted anti-tipper. For safe operation, the front anti-tipper must be installed correctly and maintained in proper condition.

⚠ CAUTION

Risk of accidents and injury due to improper installation of the Recaro® seat. During installation, ensure that the seat adapter is securely attached to the retainer bushings and that the anterior locking bolts have snapped into place so that they touch the key ring.

INFORMATION

Prior to using the power wheelchair, all the necessary mechanical adaptations (e. g. mounting special controls or a push button module) and software settings (e. g. programming the control unit) must be made to comply with the individual requirements and abilities of the user. Settings may only be made by personnel that has been authorised and trained by Ottobock.

INFORMATION

When installing special controls or a push button module, the user is obliged to have a specialist dealer give him or her detailed instructions on how to use it.

INFORMATION

The installation position of the fixture to adjust the seat height and / or seat may not be modified.

2.5 Safety requirements for operation

⚠ WARNING

Risk of accidents and injury due to defective safety functions. The user is obligated to ensure that the power wheelchair and its safety functions are in safe and proper condition before every use.

The power wheelchair may only be operated if all safety functions, e.g. the automatic brakes, are functional. Brake failure can result in serious accidents with fatal injuries.

⚠ WARNING

Risk of accidents and injury due to incorrect configuration settings. Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents.

Always test the driving characteristics of the power wheelchair after configuration / programming is complete.

Programming must only be completed by authorised personnel. Neither Ottobock nor the control unit manufacturer are liable for damages (especially in combination with special controls) caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.

⚠ WARNING

Risk of accidents and injury due to a lack of brake functionality when the brake is disengaged. Note the lack of brake functionality when the brake is disengaged, especially when moving the power wheelchair on slopes.

⚠ WARNING

Risk of injury due to tipping during operation. The C1000 DS has been approved for ascending or descending inclines up to a maximum of 17 %. Navigating inclines above this percentage value is not permitted.

The critical obstacle height of the C1000 DS is 10 cm. Crossing obstacles higher than 10 cm is not permitted. Obstacles must not be crossed while ascending or descending inclines.

Navigating stairs is not permitted.

⚠ WARNING

Risk of accidents due to tipping during operation. Reduce speed when driving downhill (e.g. select speed level 1). Obstacles must not be crossed while ascending or descending inclines. Avoid getting into or out of the wheelchair on slopes.

Before driving up slopes or over obstacles, a tilted or raised seat must be lowered to its basic position and the back be brought into the upright position. Tilting the seat back slightly is recommended when navigating downhill. Obstacles such as steps or low curbs may only be crossed at reduced speed (max. 3 km/h). Always approach obstacles at a right angle and cross over them without stopping.

⚠ WARNING

Risk of tipping when driving on unsuitable terrain. Driving on very smooth areas (e.g. icy surfaces) or extremely rough terrain (e.g. broken rock or rubble) is not permitted.

△ WARNING

Risk of tipping when using lifting platforms. Turn off the control unit of the power wheelchair when using lifting platforms, lifts, buses or trains. Engage the brake.

⚠ WARNING

Risk of tipping due to changes in the centre of gravity. The effects of changes in the centre of gravity on the behaviour of the power wheelchair on inclines, side hills or when crossing obstacles must be tested with the secure support of an attendant prior to initial use.

⚠ WARNING

Risk of accidents and injury due to incorrect lifting. Attendants may only lift the wheelchair by securely mounted components, but never by the footrests or arm rests.

⚠ WARNING

Risk of crushing or pinching in the power seat adjustment and lifting range. When activating the fixture to adjust the seat height, back angle and seat tilt, there are pinch and shear points in the area between the seat frame and wheelchair frame due to the construction. All attendants must be informed about this. In order to prevent injuries, ensure that body parts such as the hand and feet are always kept out of the danger zone, that no interfering objects such as clothing or obstacles are in the danger zone and that no unauthorised persons are present there.

⚠ WARNING

Risk of accidents and injury due to improper use of the seat height adjustment function.

■ The seat height adjustment function may only be used on a level surface.

- When using the seat height adjustment function, no unauthorised persons may be present in the danger zone.
- There must not be any interfering objects or obstacles in the adjustment range.
- All attendants must be informed that there are pinch points caused by the construction in the area between the seat frame and wheelchair frame.
- Users and attendants must not reach into the danger zone.

⚠ CAUTION

Risk of crushing or pinching in the pivot range of the footrests. Ensure limbs are clear of the danger zone when folding the footplates up or down.

⚠ CAUTION

Risk of accidents when driving without experience. Driving without experience can lead to tipping or other dangerous situations. Practice using the wheelchair on level, straightforward terrain before using it.

⚠ CAUTION

Risk of accidents due to uncontrolled driving characteristics. Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunctions. In this case, please contact your authorised dealer immediately. If any faults, defects or other hazards that can lead to personal injury are detected, the power wheelchair must be taken out of service immediately.

⚠ CAUTION

Risk of accidents and injury due to getting into or out of the wheelchair improperly. Turn off the wheelchair controls in order to get into or out of the wheelchair. Note that the arm rests are not capable of bearing full body weight, and therefore must not be used for getting into or out of the wheelchair.

⚠ CAUTION

Risk of accidents and injury due to rolling. Releasing the brake may result in uncontrolled rolling of the power wheelchair. Therefore, ensure that the brake is engaged after parking the power wheelchair.

⚠ CAUTION

Risk of accidents due to uncontrolled driving characteristics. Interference caused by electromagnetic fields may occur during the operation of the power wheelchair. The following particularities must be observed during operation:

- The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). Turn off all mobile devices when operating the wheelchair.
- The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Turn off the control unit whenever you do not need it.
- The power wheelchair has been tested according to EMC regulations.

⚠ CAUTION

Risk of accidents due to bad tyres. Visually inspect the tyres for sufficient tread depth and correct pressure before each use.

Avoid parking outside and in direct sunlight (UV radiation) since this causes the tyres to age quickly. As a result, the tread surface hardens and corner pieces break out of the

tread. Ottobock recommends replacing the tyres at least every 2 years regardless of the amount of wear.

If the wheelchair is parked for an extended period of time or the tyres overheat (e.g. in the vicinity of radiators or in case of exposure to strong sunlight behind glass), the tyres may become permanently deformed. Therefore, sufficient clearance should always be maintained from sources of heat and the C1000 DS power wheelchair should be moved frequently or placed on blocks for storage.

⚠ CAUTION

Risk of accidents due to improper clothing. Ottobock advises users to wear light-coloured clothing or clothing with reflectors during hours of darkness.

⚠ CAUTION

Risk of injury due to falling out of the wheelchair. It is strongly recommended that the user for his/her own safety always fastens the safety belt. The safety belt helps to stabilise the person sitting in the power wheelchair.

NOTICE

Risk of damage due to improper use of the seat adjustment functions. The C1000 DS power wheelchair may be equipped with seat height adjustment, back angle adjustment and seat tilt functions. The following particularities must be observed during operation:

- The seat function actuators are not intended for continuous use but only for short-term limited operation (10% use, 90% idle time). General recommendation: At maximum load capacity, 10 seconds of operation should be followed by 90 seconds of idle time. The electric seat functions are considered independently of the driving function for this purpose.
- The maximum load capacity of the seat with seat height adjustment function is 130 or 200 kg.
- The seat height adjustment function may not be operated in case of errors or malfunctions.
- If creep speed is not activated when the seat height adjustment function is operated, please visit a specialist dealer immediately. The C1000 DS power wheelchair may only be used with the seat in its lowest position until the problem is corrected.

NOTICE

Risk of damage due to excessive heat or cold. The C1000 DS can only be operated in the temperature range from -25 °C to +50 °C. It must not be operated at temperatures outside this range.

NOTICE

Risk of damage due to overloading / excessive strain. The maximum load capacity for the standard version of the C1000 DS power wheelchair is 140 kg. It is possible to have the power wheelchair retrofitted for a maximum load capacity of 200 kg. If the power wheelchair is equipped with a seat height adjustment unit, the maximum load capacity is 130 kg, optional 200 kg.

NOTICE

Damage to tyres due to excessive tyre pressure. Do not exceed the tyre pressure specified in Section 10 "Technical Data".

NOTICE

Electric interference due to electromagnetic fields. The power wheelchair has been tested according to EMC regulations. Take care of the following particularities during the operation:

The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). Please switch off all mobile devices when driving.

The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Therefore, switch off the controls whenever you do not need them.

Despite of compliance with all applicable EMC directives and standards, the power wheelchair can be affected by interference from other electric devices (e.g. department store EAS systems) or cause interference to such devices. If you notice such a behaviour, move your power wheelchair outside the interference range.

INFORMATION

Road traffic regulations must be observed when driving in road traffic.

INFORMATION

Every time you switch on the control unit, it will return to the previously selected speed level. If desired, the starting speed level can be changed in the parameter settings. These settings may only be established by personnel that has been authorised and trained by Ottobock.

INFORMATION

If the driving mode is changed while driving, the power wheelchair will accelerate or decelerate.

INFORMATION

The control unit of the power wheelchair is protected according to protection class IP 54 and the push button module according to protection class IPX4. Accordingly both can be used in bad weather (e. g. rain). The control unit or push button module is approved for indoor and outdoor operation and meets the requirements regarding environment and splashed water.

INFORMATION

The control unit of the power wheelchair must be reactivated after an emergency stop. In the event of communication problems in the bus system of the control unit, the system triggers an emergency stop and thus prevents any uncontrolled functions. If the driving function is still not available after reactivating the controls, the brake must be released to switch over to pushing mode. In this case, contact a specialist dealer as soon as possible.

INFORMATION

If the control unit does not emit an error signal as soon as the brake lever is activated, this indicates a malfunction. The settings must be checked by a specialist dealer.

INFORMATION

Manoeuvring is only allowed at reduced speeds.

INFORMATION

During the use of the power wheelchair, electrical discharges may occur which are caused by factors such as friction (high voltage with low current; discharge via the user). However, these do not represent a health hazard.

This reduction in comfort – under consideration of the environmental conditions the user is exposed to – can be eliminated by custom fabrication measures (installing a mechanical discharge contact / ground strap to the wheelchair frame).

Electrostatic discharge may also occur if the power wheelchair is equipped with puncture-proof tyres. Retrofitting the wheelchair with pneumatic tyres can correct this problem.

2.6 Safety requirements for care, maintenance and disposal

⚠ CAUTION

Risk of accidents and injury due to improper maintenance, repairs or adjustments. Maintenance work on the power wheelchair may only be completed by personnel authorised and trained by Ottobock. This also applies to all repairs and settings on the brake. Incorrect settings can lead to a loss of braking power.

⚠ WARNING

Risk of injury due to explosive gases. Explosive gases can develop while the batteries are charging. The following safety instructions must be observed under all circumstances:

- Ensure sufficient ventilation when charging the batteries in an enclosed room.
- Smoking and fire are not permitted.
- Sparks must be avoided.
- Do not cover the air vents in the trim.

⚠ CAUTION

Risk of injury due to uncontrolled movements. The fuse must always be removed for any maintenance work where the battery cover is open.

NOTICE

Unauthorised battery replacement. Battery replacement or changing the battery installation location may only be completed by a specialist dealer. The charging profile of the battery charger established at the factory corresponds to the battery included in the scope of delivery and may not be independently altered. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

NOTICE

Risk of damage to the electronics caused by short circuits. Always use insulated tools during all maintenance work on the batteries. Ensure correct polarity is used when connecting the battery cables. The black cable must be connected to the negative pole and the red cable to the positive pole.

NOTICE

Risk of damage to electronics due to water penetration.

Never use a water jet or high-pressure cleaning apparatus to clean the power wheelchair. Water must not come into direct contact with the electronics, motor or battery.

NOTICE

Risk of damage due to uncontrolled movements. When raising the wheelchair, secure it against slipping and tipping sideways by using a suitable block under the drive unit sustainer.

INFORMATION

The power wheelchair must be inspected for functional reliability and operational safety at least once a year by an authorised specialist.

INFORMATION

When closing the battery drawer, the release bolt must engage in the recess provided for this purpose.

INFORMATION

Piston rods do not require lubrication. They are maintenance-free.

INFORMATION

Defective batteries must be disposed of properly in accordance with country-specific regulations.

2.7 Requirements for the user

⚠ CAUTION

Risk of accidents and injury due to improper operation.

The power wheelchair may only be operated by properly informed users. For this purpose, the user and any attendants must receive instruction in the use of the power wheelchair from specialists trained and authorised by Ottobock. The operator must have read and understood all the information in the instructions for use. The wheelchair may not be operated in cases of exhaustion or under the influence of alcohol or medications.

The user must not have any mental limitations which can temporarily or permanently restrict attentiveness and judgment in road traffic.

2.8 Safety functions

INFORMATION

In dangerous situations, the C1000 DS can be turned off at any time using the on/off button. When the button is pressed, the power wheelchair brakes immediately and the electric functions will be stopped.

Malfunctions such as an insufficient supply of power to the brake are recognised by the software, triggering the emergency stop brake or reducing the speed of the power wheelchair. At the same time, a warning signal will be emitted.

INFORMATION


The control unit of the C1000 DS power wheelchair must be reactivated after an emergency stop. In the event of communication problems in the bus system of the control unit, the system triggers an emergency stop and thus prevents any uncontrolled functions. If the driving function is

still not available after reactivating the controls, the brake must be released to switch over to pushing mode. In this case, contact a specialist dealer as soon as possible.

INFORMATION

The control unit of the power wheelchair switches to a safe mode at elevated temperatures and after driving uphill for extended periods of time, limiting the performance of the power wheelchair. However, the user is able to drive the power wheelchair out of a hazardous situation at any time. Once the excessive temperature has been alleviated (this may take several minutes depending on the ambient temperature), the power wheelchair is once again fully functional.

2.9 Warning symbols and type plates

| Label/plate | Meaning |
|-------------|--|
| | Risk of pinching. Do not reach into the danger zone. |

Table 1 Signage on the C1000 DS

3 Product description

The C1000 DS power wheelchair is suitable for indoor and outdoor use. The powerful drive system with two 12 V gel batteries combined with spring-mounted drive wheels allows obstacles to be crossed easily (Category B according to EN 12184) and offers safe operating performance.

The power wheelchair is controlled by an enAble50 wheelchair control system. The latter includes a control panel to enter driving commands and display the current status as well as a controller that controls the drive motors and other electrical functions based on the input data. Data transmission is realised via a CAN-bus system.

The enAble50 is programmable, which means the control unit can be customised to individual user requirements. For example the speed, acceleration and delay values can be adjusted.

The two steering casters are controlled by separate steering motors. This considerably improves the driving characteristics, especially on uneven terrain or at higher speeds. The steering casters driven by steering motors turn by 90° or 360° depending on the programming. This results in excellent manoeuvrability, both indoors and outdoors.

Special features of the C1000 DS include:

- Serviceability due to easy, straightforward access to all components.
- Individual adaptation possibilities with options and custom fabrication using modular components (chassis, seating system, control unit, accessories).

The modular design of the C1000 DS makes it possible to add other modules and devices in addition to the main components (see Figure 1), e. g. electric seat adjustment, special controls, tray, etc. These options are described in Section 7.

4 Delivery and preparation for use

4.1 Delivery

INFORMATION

The options included in the scope of delivery depend on the product configuration that has been selected for the power wheelchair.

The scope of delivery includes:

- Pre-set power wheelchair with main components (see Figure 1)
- Charger
- Instructions for use
- Options (see Section 7)

The specialist dealer delivers the power wheelchair ready for use. All settings correspond to the indications on the order form or are adjusted directly on site by the dealer. The power wheelchair is adapted to the personal requirements of the individual.

The functionality of the individual components can be tested by following the instructions in Section 6. Section 8 includes a description of possible malfunctions.


Figure 1 Main components

- 1 Back rest
- 2 Joystick and control unit
- Arm rest (side panel)
- 4 Seat cushion

- 5 Wheel lock release
- 6 Footrest
- 7 Motor with drive wheel


Figure 2 Anti-tipper

1 Anti-tipper front, both sides

4.2 Settings

⚠ WARNING

Risk of accidents and injury due to incorrect settings. All subsequent adaptation or adjustment work may only be carried out by authorised specialist dealers.

The following tools are required to adapt the wheelchair to the needs of the user and / or to complete the maintenance tasks:

- Allen wrenches, sizes 3 to 6
- Ring or open-end wrenches, sizes 8, 10,13
- Ratchet and socket, size 17
- Torque wrench
- Phillips screwdriver
- Screwdriver for pan head screws
- Tyre repair accessories

4.3 Initial operation

▲ WARNING

Risk of suffocation. Packaging materials must be kept out of reach of children.

All components must be checked for completeness and functionality prior to initial operation. Before the power wheelchair can be switched on, the 80 A fuse must be inserted into the fuse holder to the left rear of the battery case (see Figure 3).

To do so, remove the fuse from the protective cover on the control panel, open the cover of the fuse holder and insert the fuse.

Make sure that the fuse is pressed into the centre of the spring contacts and that it is not positioned at an angle. Close the cap until you feel it snap back into place.


Figure 3 Fuse

- 1 Cap of fuse housing, open
- 2 Fuse

Initial operation with electrical back angle adjustment

⚠ CAUTION

Risk of accidents and injury due to an incorrectly installed back rest. Ensure that the locking mechanism engages in the cross bolt. Verify that the cross bolt and locking mechanism are securely in place.

Folding up the back rest:

- 1. Fold the back rest up.
- 2. Insert the cross bolt at the end of the drive unit into the holder.
- 3. Push the lever down. The locking system engages in the cross bolt.

Figure 4 Inserting the cross bolt into the holder; locked cross bolt

- 1 Cross bolt
- 2 Locking system with lever

5 Transportation and storage

5.1 Transportation in a Wheelchair Accessible Vehicle

⚠ WARNING

Risk of accidents and injury due to improper use for transportation in a wheelchair accessible vehicle. Use in a wheelchair accessible vehicle is not permissible if the following options are used:

- Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system nn
- 94 Ah (C5) gel batteries nn
- The maximum load capacity is limited to 100 kg.

⚠ WARNING

Risk of accident and injury due to improper use for transportation in a wheelchair accessible vehicle. The power wheelchair may only be used for transportation in a wheelchair accessible vehicle if the safety components (such as lap belt) which are offered by Ottobock and suitable restraint systems are used.

Only one person may be transported in the C1000 DS at any one time. For optimum protection of the occupants in

case of an accident, use the seats and corresponding restraint systems installed in a wheelchair accessible vehicle while it is in operation.

For more information on the use of the C1000 DS as a seat in a wheelchair accessible vehicle, please refer to our brochure "Using your Wheelchair / Mobility Base with Seating Shell or Buggy for Transportation in Wheelchair Accessible Vehicles", order number 646D158.

Ensure that the wheelchair is fully secured with tensioning straps during transportation in a wheelchair accessible vehicle.

The frame has four fixing eyelets for the attachment of straps (see Fig. 5).

Before transporting the power wheelchair, switch off the control unit and engage the wheel lock.

5.1.1 Required Accessories

To use the power wheelchair as a seat in a wheelchair accessible vehicle, additional accessories have to be mounted (491S00=SK028 Anchor Point Kit).

The qualified personnel which fit the wheelchair can provide more information.

5.2 Transport of the Power Wheelchair

⚠ WARNING

Risk of accidents and injury due to loosening plug and screw connections. Screw and plug connections must be loosened for all adjustment work and when folding the back rest for transportation. This can cause uncontrolled movements of wheelchair components. To avoid injuries at pinch and shear points, carry out such tasks with the support of another person. Ensure that body parts such as the hands or head are kept out of the danger zone at all times.

⚠ CAUTION

Risk of injury due to insufficient restraints. Secure the C1000 DS power wheelchair during transportation in another vehicle by tightening the tensioning straps properly. Only attach tensioning straps to the corresponding eyebolts and specified mounting points.


⚠ CAUTION

Risk of pinching or crushing in the pivot range of the footplate(s). Ensure limbs are clear of the danger zone when folding the footrests up or down.

NOTICE

Risk of damage due to falling. The maximum net weight of the C1000 DS power wheelchair is 129 kg. Hoisting devices and transportation vehicles used for transportation must have sufficient capacity.

The eyebolts at the font and back of the power wheelchair serve as mounting points for tensioning straps.


Figure 5 Eyebolts front, eyebolts rear

- 1 Eyebolts, front
- 2 Eyebolts, rear

5.2.1 Preparing for transportation

For the purpose of transportation, the size of the wheelchair can be reduced by folding the footplates up or removing them, removing the side panels and folding down the back rest.

Before transporting the C1000 DS power wheelchair, switch off the control unit and engage the wheel lock.

Folding the back rest

1. Remove the side panels (see Section 6.1.3).

- Pull the release strap to the rear.
- With your hand, fold the back rest down onto the seat bottom.

Figure 6 Folding the back rest

Folding the back rest with electric back angle adjustment

- 1. Remove the side panels (see Section 6.1.3).
- 2. Unlock the cross bolts on the lower end of the drive unit by pushing the lever up (see Section 4.3 "Initial operation", Figure 4).
- 3. Remove the cross bolt from the holder.
- With your hand, fold the back rest down onto the seat bottom.

Figure 7 Folding size

Prior to operating the power wheelchair with back angle adjustment, proceed in the reverse sequence.

5.3 Storage

The C1000 DS power wheelchair must be stored in a dry place.

For transportation and storage, an ambient temperature between -40 $^{\circ}$ C and +65 $^{\circ}$ C must be observed.

INFORMATION

If your power wheelchair is not moved for several days, permanent colour changes may occur where the wheelchair comes into contact with the surface it is standing on. A suitable mat should be used when parking the wheelchair for extended periods of time.

INFORMATION

Tyres contain chemical substances that can react to other chemical substances (such as cleaning agents, acids, etc.). Black tyres contain soot particles. They can leave black traces on the ground. For this reason, Ottobock recommend choosing grey tyres when using the power wheelchair mainly indoors.

INFORMATION

Direct exposure to sunlight / UV radiation causes the tyres to age more quickly. As a result, the tread surface hardens and corner pieces break out of the tread. Avoid parking the wheelchair outdoors whenever possible

INFORMATION

The tyres should be replaced every 2 years regardless of wear and tear.

INFORMATION

Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

6 Operation

6.1 Adjustment possibilities

⚠ WARNING

Risk of accidents and injury due to incorrect settings. All subsequent adaptation or adjustment work may only be carried out by authorised specialist dealers.

△ WARNING

Risk of accidents and injury due to loosening plug and screw connections. Screw and plug connections must be loosened for all adjustment work and when folding the back rest for transportation. This can cause uncontrolled movements of wheelchair components. To avoid injuries at pinch and shear points, carry out such tasks with the support of another person. Ensure that body parts such as the hands or head are kept out of the danger zone at all times.

⚠ CAUTION

Risk of accidents due to unsecured screw connections.

When screw connections with thread locks are loosened, they must be replaced with new ones or secured with a medium-strength liquid thread lock substance (e.g. Loctite 241°).

After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened. Observe the applicable torque specifications.

Various adjustments can be made to the power wheelchair.

The seat height, seat width and seat angle have been set in accordance with the information on the customer order form and may only be changed by the specialist dealer.

The standard springs in the drive wheels and the springs of the steering caster forks must be adjusted by a specialist dealer.

The following items can be adapted by the user:

- Back angle
- Height and position of the arm rests
- Lower leg length

The footrests and side panels can be removed if needed.

6.1.1 Back rest

The back rest can be set to 4 different positions by pulling the release strap (base of the back rest, Figure 8).

Once the back rest is in the desired position, release the strap to engage the locking mechanism.

Figure 8 Release strap for the back rest

6.1.2 Adjusting the luggage carrier

The luggage carrier may be adjusted to one of four different positions.

1. Slightly raise the back of the luggage carrier.

Pull or push the luggage carrier until it engages in the next notch.

Figure 9 Adjusting the luggage carrier

6.1.3 Side panel with arm rest

Removing the side panel with arm rest

- 1. Loosen the thumb screw at the lower end of the side panel holder.
- 2. Pull the side panel with arm rest up and out.
- 3. After remounting the side panel, do not forget to firmly retighten the thumb screw.

If the control panel is attached to the arm rest, it must be removed before taking out the side panel (Section 6.1.4).

Figure 10 Loosening the thumb screw for removal of the side panel with arm rest

Adapting the arm rests to the upper arm length

- 1. Loosen the Allen head screw using an Allen wrench size 3.
- Move the side panel with arm rest up or down along the rail into the desired position.
- 3. Retighten the Allen head screw.

Figure 11 Adapting the arm rest to the upper arm length

Adapting the arm rests to the forearm length

- Loosen the two Allen head screws on the underside of the arm rest using an Allen wrench size 3. Three set screws must be loosened on the side where the control panel is attached (see Figure 12).
- 2. Move the arm rest with the control panel holder or replacement tube forward or backwards to the desired position along the attachment rail.
- 3. Retighten all set screws.

Figure 12 Adapting the arm rests to the forearm length

6.1.4 Control panel

⚠ CAUTION

Risk of accidents and injury due to uncontrolled movements of the wheelchair. Before changing the position of the control panel or installing/ removing it, turn off the control unit of the power wheelchair.

NOTICE

Cable damage. Positioning the cables incorrectly can lead to pinching and damage to the cables. The cables must not be attached too tightly or too loosely. Avoid bending or squeezing the cables.

Adapting the control panel to the arm length

The control panel is attached to a rail which is located beneath the arm rest.

- Loosen the three set screws on the underside of the arm rest using an Allen wrench size 3.
- Move the arm rest with the control panel holder or replacement tube forward or backwards to the desired position along the attachment rail.
- Retighten all set screws.

Figure 13 Adapting the control panel

INFORMATION

If the control panel rail is too long, the protruding part can simply be cut off with a saw.

Changing the control panel mounting side

The standard configuration has the control panel mounted on the right side. If desired by the user, it can also be mounted on the other side of the power wheelchair.

INFORMATION

Only authorised personnel may change the control panel mounting side.

6.1.5 Footrests

⚠ CAUTION

Risk of pinching. Ensure limbs are clear of the danger zone when folding the footrests up or down.

NOTICE

Risk of damage to the front lighting. Do not swing the footrests to the outside or damage to the front lighting may occur.

INFORMATION

Complete removal of the footrest from the adapter is only permitted for transportation of the C1000 DS power wheelchair.

Removing the footrests

- Remove the calf band.
- Fold the footrests up, release the footrest locking mechanism (see Figure 14) and swing the footrest to the inside.
- 3. Pull the footrest up to remove it.

Installing the footrests

- Reinsert the footrest from above by placing it in the adapter and swinging it forward until the footrest locking mechanism engages.
- 2. Reattach the calf band to the adapter.

Figure 14 Footrest locking mechanism

Adapting the footrests to the lower leg length

- 1. Loosen the screws on the footplate bar using an Allen wrench size 4 (see Figure 15).
- Move the footrests up or down to adapt the height to the individual lower leg length and seat cushion thickness.

3. Retighten the screws.

INFORMATION

When making the adjustment, ensure that the footplate bar (Figure 15, item b) slides at least 60 mm into the swivel segment (Figure 15, item c). After completing the adjustment, tighten the set screw with a torque of 8 Nm.


Figure 15 Adjusting the lower leg length

6.2 Getting into and out of the power wheelchair

⚠ CAUTION

Risk of injury if the power wheelchair starts rolling. Switch off the control unit before getting into and out of the power wheelchair. This will automatically engage the motor brake.

NOTICE

Risk of breaking the arm rests due to excessive strain. Note that the arm rests are not capable of bearing full body weight, and therefore must not be used for getting into or out of the wheelchair.

The modular design of the power wheelchair and the ease with which you can remove the side panels and footrests make it easy to get into and out of the wheelchair from the side or from the front.

Users can get into and out of the wheelchair using the best method for each individual.

6.2.1 From/to the side

Getting into the wheelchair from the side requires the user to remove the right or left side panel depending on the side where

the user will get into the wheelchair. Bring the power wheelchair as close as possible to the place where its user is sitting.

If the control panel is located on the side where the user is getting in or out, carefully remove the control panel. Remove the side panel as described in Section 6.1.4 (remove the footrest if necessary; see Section 6.1.5).

The user can then slide onto the wheelchair seat from the side. Use of a transfer board will make this even easier.

NOTICE

Risk of damage to the front lighting. Do not swing the footrests to the outside or damage to the front lighting may occur.

6.2.2 From / to the front

Flipping up the two footplates makes it possible for the user to get into and out of the wheelchair from the front (Figure 16). If required, it is also possible to remove the footrests from their adapter after swinging them out of the way (Figure 17; please observe the safety information in Section 6.1.5).

The assistance of an attendant or a transfer lifter makes it easy for the user to get into or out of the power wheelchair. A rotation plate offers support.

Figure 16 Footrests folded up

Figure 17 Footrests removed

6.3 Control unit

Possible impairment of the power wheelchair driving characteristics. The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). For this reason, please switch off all mobile devices when driving.

Risk of damage to other devices. The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Turn off the control unit whenever you do not need it.

6.3.1 Control panel

The power wheelchair is operated using the control panel.

The control panel consists of a keypad, LCD monitor and joystick. The programming receptacle and two inputs for external push-buttons are located on the underside. The power wheelchair is turned on and off using the keypad; driving commands can be entered and the status of certain functions and components can be displayed.


Figure 18 Control panel

- Joystick
- 2 Direction indicator lights, left
- 3 Warning flasher
- 4 LCD monitor
- 5 On/Off button

- 6 Mode button
- 7 Horn
- B Direction indicator lights, right
- 9 Light


Figure 19 Control panel connections on the underside

- 1 External On/Off button connector
- 2 Connector for external mode button
- 3 Programming receptacle

On/Off button

The On/Off button is used to turn the power wheelchair on and off.

Mode button

Briefly pressing the mode button increases the speed level. After reaching the maximum speed level, you can change back to speed level 1 by pressing the button again (1-2-3-4-5-1-...). Pushing and holding the button for at least 2 seconds

switches the control unit to the "Electric seat function" menu.

Horn

The horn will sound as long as the horn button is being pressed.

Warning flasher

Pressing the warning flasher button activates or deactivates all four direction indicator lights.

Light

The front and rear lights are activated and deactivated by pressing the light button.

Direction indicator lights

The right or left front and rear direction indicator lights are activated and deactivated by pressing the "direction indicator light, right" or "direction indicator light, left" buttons respectively. The direction indicator lights automatically turn off after 20 seconds.

LCD monitor

The LCD monitor is the communication interface between the user and the control unit. It indicates the selected speed level, the remaining battery capacity and the status of electric op-

tions and special functions as well as warnings and errors. All display symbols are visible during the initialisation phase.


Figure 20 LCD monitor with all symbols

- 1 Direction indicator lights, left
- 2 Battery capacity
- 3 Electric back rest
- 4 Electric seat tilt
- 5 Control panel
- 6 Drive-away lock
- 7 Speed level
- 8 Open-end wrench

- 9 Excess temperature
- 10 Direction indicator lights, right
- 11 Wheel lock for drive wheels
- 12 Warning
- 13 Creep speed
- 14 Power module
- 15 Drive motor
- 16 Light

6.3.2 Separate LCD monitor

A separate LCD monitor with infrared and Bluetooth (optional) is available to operate special controls and for environment control functions. Please see Section 7.9.1 for more information.

6.3.3 Switching on and off

▲ WARNING

Danger to life due to lack of brake functionality. Please note that the wheel lock release lever must be locked prior to operating the power wheelchair. The automatic wheel lock must be operational and functional.

⚠ CAUTION

Risk of accidents due to incorrect tyre pressure. Visually inspect the tyres of the power wheelchair for sufficient tread depth and sufficient tyre pressure before each use. Incorrect tire pressure reduces the life of the tyres and impairs the driving characteristics.

The control unit of the power wheelchair is switched on and off by pressing the On/Off button (see Figure 18, item 5). If the control unit is not used for a certain period of time, which is adjustable, the wheelchair is switched off automatically. It is also possible to switch off the power wheelchair with the On/Off button while driving. In this case, the wheelchair brakes immediately until it stops.

INFORMATION

Under the parameter settings, the specialist dealer can establish the default speed level or menu for the power wheelchair after it is turned on according to the wishes of the user.

6.3.4 Driving function

⚠ CAUTION

Risk of injury due to uncontrolled operation of the power wheelchair. Switch off the control unit on the power wheelchair whenever you do not need the driving function. This will prevent you from unintentionally activating the joystick.

⚠ CAUTION

Risk of injury due to tipping during operation. Observe the following safety information when operating the power wheelchair:

- Operation of the wheelchair is limited to inclines up to a maximum of 17 %.
- When driving downhill, reduce speed in accordance with the degree of the incline.
- Do not cross obstacles higher than 10 cm.
- Do not cross steps or curbs without reducing speed.

INFORMATION

The control unit of the power wheelchair switches to a safe mode at elevated temperatures and after driving uphill for extended periods of time, limiting the performance of the power wheelchair. However, the user is able to drive the power wheelchair out of a hazardous situation at any time. Once the excessive temperature has been alleviated (this may take several minutes depending on the ambient temperature), the power wheelchair is once again fully functional.

The C1000 DS power wheelchair has been approved for ascending or descending inclines of up to a maximum of 17%. Driving on steeper inclines or slopes is not permitted.

In order to navigate downhill grades safely, the speed must be reduced according to the slope (e.g. select speed level 1).

The critical obstacle height of the power wheelchair is 10 cm. Crossing obstacles higher than 10 cm is not permitted. Obstacles such as steps or low curbs must be crossed at reduced speeds.

If there are obstacles in the travel path, it is important to drive around them with a large margin of safety. On uneven ground, the driving characteristics of the wheelchair may get out of control. Therefore the speed must always be adapted to the ground conditions.

The battery indicator and speed level are shown in the driving menu on the LCD monitor.

Use the joystick for driving. The further the joystick is moved away from the centre position, the faster the power wheel-chair will drive in this direction.

The maximum speed with full deflection of the joystick depends on the selected speed level. Releasing the joystick automatically activates the brake function, which brings the

wheelchair to a halt. When standing still, the mechanical wheel locks are automatically active so the power wheelchair cannot roll.

The power wheelchair has a programmable number of speed levels (delivery condition = 5 speed levels). Use the "M" button to increase the speed level. After reaching the maximum speed level, the control unit switches back to speed level 1.

Adapting the operating performance

⚠ WARNING

Risk of accidents and injury due to incorrect configuration settings. Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents. Always test the driving characteristics of the power wheelchair after configuration / programming is complete. Programming must only be completed by authorised personnel. Neither Ottobock nor the control unit manufacturer are liable for damages (especially in combination with special controls) caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.

The speed, acceleration and deceleration values can only be adjusted to individual user requirements by the specialist dealer via a handheld programming device that is connected to the control panel.

6.3.5 "Battery Capacity" indicator

INFORMATION

Immediately after switching the power wheelchair on, the battery indicator shows the battery capacity that was saved before the wheelchair was switched off the last time. The exact battery status is displayed after approximately 2 minutes.

The battery indicator on the LCD screen is divided into 7 segments that show the current battery capacity.

A charge of 100% corresponds to 7 segments on the battery symbol. As the remaining battery capacity decreases, the LCD segments turn off one by one (see Table 2).

If all the segments have disappeared and the battery symbol is flashing, the battery is under voltage. Since any further use will lead to battery damage, a warning signal also appears.

If the battery symbol and all the segments are flashing, this means that the battery is over voltage. Since this will lead to battery damage, a warning symbol also appears.

The charging process is indicated by the battery segments lighting up one after the other. When the battery is charging, the driving function is locked out.

| Display symbol(s) | Information |
|-------------------|--|
| | Driving menu with speed level and battery capacity |
| | Low battery capacity |
| | Charging process with drive stop |
| | Battery is under voltage, with warning symbol |
| △ | Battery is over voltage, with warning symbol |

Table 2 Battery indicator on the control panel

6.3.6 Controlling additional electrical functions

INFORMATION

For more information on the additional electrical functions, please see Section 7.

Additional electrical functions, e.g. electric back angle adjustment and seat tilt, are accessed by pressing the "M" button (approximately 2 seconds).

To change between the different functions, briefly press the button "M" or move the joystick to the right. The function that is currently selected is shown on the LCD monitor. The respective function can be extended or retracted by moving the joystick ahead or back. The electric motor moves the respective function as long as the joystick is deflected and stops at the end positions.

Electric back angle adjustment: When the joystick is deflected ahead or back, the back rest moves ahead or back accordingly. The electric motor moves the back rest as long as the joystick is deflected and stops at the end positions.

Electric seat tilt: Deflecting the joystick to the rear will electrically tilt the seat to the rear. Deflecting the joystick to the front tilts the seat ahead until it is level (factory setting, may

be adapted by the specialist dealer).

Seat height adjustment fixture: Deflecting the joystick ahead or back raises or lowers the electric lift seat. As soon as the seat height adjustment fixture leaves the lowest position, the speed is reduced (creep speed) because the power wheelchair is less tip-resistant. On the control panel, this is indicated by the snail symbol on the LCD monitor (see Figure 20).

Coupled seat adjustment: Deflecting the joystick ahead or back moves the back rest and seat forward and down or back and up simultaneously.

Electric footrests: Deflecting the joystick ahead or back raises or lowers the electric footrests.

Pressing the "M" button (approximately 2 seconds) switches back to driving mode.

| Display symbol(s) | Information |
|-------------------|--------------------------------|
| 5-5 | Electric back angle adjustment |
| 5-5 | Electric seat tilt |

| Display symbol(s) | Information |
|-------------------|--|
| | Electric seat height adjustment |
| | Coupled electric seat adjustment (back rest and seat tilt) |
| | Electric footrests, coupled |
| | Electric footrest, left |
| | Electric footrest, right |
| | Creep speed |
| | Attendant control |

Table 3 Display of electric seat functions on the control panel

6.3.7 Drive-away lock

INFORMATION

The drive-away lock can only be activated in the manner described below if the "Drive-away lock" parameter is set to "On". At the factory, the drive-away lock parameter is set to "Off". The factory setting may be set to "On" by the specialist dealer or at the factory according to the order. Please ask your specialist dealer about the setting selected for your wheelchair.

The control unit of the C1000 DS power wheelchair features an electric drive-away lock to disable the driving function.

Activation via control panel

- 1. With the control unit turned on, press and hold the mode button for at least 5 seconds.
- A short beep confirms activation of the drive-away lock (if the feature that returns a beep with each push of a button is activated, then the second beep confirms activation of the drive-away lock).
- 3. The control unit is turned off automatically.

Activation of the drive-away lock is indicated by the key symbol on the display.

| Display symbol(s) | Information |
|-------------------|-----------------|
| • | Drive-away lock |

Table 4 Display of the drive-away lock on the control panel

Activation via separate LCD monitor

When using one of the special controls, the drive-away lock is activated via the separate LCD monitor. Select the menu item "Drive-away lock" in the menu "Settings". Activation is confirmed with an acoustic signal and the control unit is switched off.

Figure 21 Information window with drive-away lock

Deactivation via control panel

- 1. Press the On/Off button to turn on the control unit.
- 2. Move the joystick to the front end position until you hear a beep.
- 3. Move the joystick to the rear end position until you hear a beep.
- 4. Release the joystick.
- Another beep confirms the successful deactivation of the drive-away lock.
- The speed level and battery indicator appear on the LCD monitor and the control unit is ready for driving.

Now the drive-away lock is deactivated and the wheelchair can be driven.

INFORMATION

If the joystick is not moved correctly, the drive-away lock remains activated. In order to deactivate the drive-away lock, the control unit must be switched off. The power wheelchair can be switched on again and the drive-away lock can be deactivated.

Deactivation via separate LCD monitor

An information window is shown on the LCD monitor after the power wheelchair is switched on. The deactivation process is the same as described under "Deactivation via control panel".

6.3.8 Steering Lock

INFORMATION

The steering lock can be activated in the manner described below if the "Steering Lock" parameter is set to "On". The steering lock parameter is set to "Off" by default. The factory setting may also have been set to "On" by the specialist dealer or at the factory according to the order. Please ask your specialist dealer about the setting selected for your wheelchair.

The enAble50 control unit in the power wheelchair may be programmed with the optional "Steering Lock" parameter.

If the steering lock is activated, the caster wheels or the drive wheels (depending on the type) are set straight ahead. This makes it easier to negotiate narrow ramps and to use elevators.

When the joystick is deflected forward, the power wheelchair drives straight ahead.

- When the joystick is deflected backward, the power wheelchair drives straight back.
- Turning is not possible when the steering lock is activated.
- The warning flasher lights flash when the steering lock is activated (standard programming).

Activating the Steering Lock

- 1. Bring the power wheelchair to a full stop.
- 2. Push the [Warning Flasher] button. The steering lock is activated. All warning flasher lights flash (standard programming).
- 3. Move the joystick forward/backward. The power wheelchair only drives straight ahead or back.

Deactivating the Steering Lock

- 1. Bring the power wheelchair to a full stop.
- 2. Push the [Warning Flasher] button. The steering lock is deactivated. All warning flasher lights stop flashing.
- 3. Move the joystick. The power wheelchair now drives in the chosen direction without steering lock.

6.4 Lighting

The C1000 DS power wheelchair comes with effective LED front and rear lighting as standard equipment. The lighting functions are activated and deactivated via the controller. The warning flashers and right / left direction indicator lights are activated using the control panel.

Front lights

INFORMATION

The front lights may only be installed and removed by authorised personnel.

An LED headlight and LED direction indicator light are attached to the right and left splash guards of the drive wheels, respectively.

Figure 22 Front lights

Rear lights

Two LED rear lights with integrated LED direction indicator lights are installed in the steering motor mouldings at the rear of the power wheelchair.

Reflectors

In addition to the rear lights, a reflector is installed on the right and left underneath the luggage carrier.

INFORMATION

When loading the luggage carrier, ensure your luggage does not cover the reflectors.

Figure 23 Rear lights, reflectors

- 1 LED rear light
- 2 LED direction indicator light
- 3 Reflectors

6.5 Releasing and engaging the wheel lock

⚠ WARNING

Danger to life due to brake malfunction. Incorrect brake settings can lead to a loss of brake functionality and therefore to serious bodily injuries and even death. Repairs and adjustments to the brake may only be carried out by authorised personnel.

⚠ WARNING

Risk of accidents and injury due to lack of brake functionality.

- No brake functionality is available when the brake is unlocked (push mode).
- When moving the power wheelchair on an incline, the person pushing must provide the appropriate brake force.
- The brake function may only be released in the presence of an attendant.
- Should the user be unable to release the brake himself, the brake can be released by the attendant.

NOTICE

Risk of damage when parking the wheelchair with the brake deactivated. Deactivating the brake may result in uncontrolled rolling of the power wheelchair. Therefore, make sure that the brake is engaged when parking the power wheelchair.

In case of control unit failure or insufficient battery capacity, the power wheelchair can be pushed. To do so, the brake is released by mechanically unlocking it. The brake release mechanism is located between the frame and seat (see Figure 24).

Releasing the brake

For safety reasons, you first have to pull out the brake release bolt.

INFORMATION

In this position, the brake still has to be engaged! If the pushing function is already active, please contact your specialist dealer.

INFORMATION

When the joystick is activated while the brake is deactivated, the control unit emits an error signal on the control panel. If this is not the case, there is a malfunction that must be corrected immediately by a specialist dealer.

Now the brake release lever has been released and can be pushed back to the stop.

Figure 24 Releasing the wheel lock

- 1 Brake release bolt
- 2 Brake release lever

In this position, the control unit will recognise that the brake has been released and automatically deactivate the driving function. As soon as the joystick is deflected, a warning will appear on the LCD monitor.

INFORMATION

Once the brake release lever has been released, all braking systems are deactivated.

INFORMATION

The control unit may be switched off in push mode.

Activating the brake

To activate the brake, push the brake release lever forward until the brake release bolt engages. To reactivate the driving function, switch the wheelchair control unit off and then back on.

6.6 Batteries

⚠ CAUTION

Using an incorrect battery type. The batteries and charger constitute a harmonised system in order to guarantee optimised charging characteristics. Using an incorrect battery type can result in health hazards due to leaking battery acid, permanent damage or destruction of the battery.

- Battery replacement or changing the battery installation position may only be completed by a specialist dealer.
- Replacement batteries must always be of the same type (Ah & V).
- The charging profile of the battery charger established at the factory corresponds to the battery included in the scope of delivery and may not be altered independently.

⚠ CAUTION

Risk of injury due to electric current. Always remove the fuse before doing any work on the batteries.

⚠ CAUTION

Risk of injury during tasks that involve lifting. The weight of a battery is approximately 20 kg. Use hoisting devices of sufficient capacity to remove or install the batteries or carry out these tasks with the aid of an assistant.

NOTICE

Risk of battery damage due to falling. Hoisting devices used for transportation must have a sufficient capacity. Secure the batteries against falling.

NOTICE

Risk of battery damage. Ensure correct polarity is used when connecting the battery cables. The black cable must be connected to the negative pole and the red cable to the positive pole.

NOTICE

Risk of property damage due to improper maintenance work on the batteries. All maintenance work may only be completed by authorised personnel.

The standard version of the C1000 DS power wheelchair is equipped with 2 maintenance-free 12 V gel batteries with a capacity of 63 Ah (C5). Higher-performance batteries with a capacity of 94 Ah (C5) are optionally available. The batteries are located beneath the seat of the C1000 DS power wheelchair.

6.6.1 Charging the batteries

⚠ WARNING

Risk of injury due to explosive gases. Explosive gases can develop while the batteries are charging. The following safety instructions must be observed under all circumstances:

- Ensure sufficient ventilation when charging the batteries in an enclosed room.
- Smoking and fire are not permitted.
- Sparks must be avoided.
- Do not cover the air vents in the trim...

NOTICE

Battery damage due to deep discharge. Extended operation with low battery capacity results in deep discharge of the battery, causing damage. When driving with the battery discharged, the control unit of the power wheelchair switches to power-saving creep speed.

The remaining battery capacity determines the range of the C1000 DS power wheelchair. The battery capacity is influenced by many factors. Besides the temperature, the age of the batteries and the time and extent of driving, the way you charge the batteries also has a considerable effect on their capacity and thus on the range.

For an optimum charging frequency, the following applies:

- The batteries can be charged at any time, regardless of the remaining battery capacity.
- It takes about 8 hours until a discharged battery (only one flashing segment) is fully charged. When the charging process is complete, the battery charger can remain connected to the power wheelchair with no risk of overcharging or damaging the battery. The battery charger features a programmed recharging phase that will maintain the battery capacity at the level that has been reached.
- For daily use, connecting the charger overnight is recommended so that the full capacity is available every day.
- If the wheelchair is not used for extended periods of time, the batteries gradually discharge. If the power wheelchair is not driven for extended periods of time, a charging cycle must be carried out at least 1 x per week in order to maintain the battery capacity. Removing the fuse is recommended for extended periods of disuse.
- Total discharge of the batteries (deep discharge) should be avoided.
- When the batteries are charging, the power wheelchair control unit must be switched off to allow all of the charging current to be fed into the battery.

The following must be observed when charging the batteries:

- Only the battery charger provided by Ottobock may be used for charging. Failure to comply with this requirement will render the warranty null and void.
- The voltage settings on the battery charger must correspond to the voltage used in your country.

6.6.2 Charger

⚠ WARNING

Explosion hazard due to sparks. Always switch off the battery charger and remove the mains plug before disconnecting the battery.

NOTICE

Risk of damage to the battery charger / risk of damage caused by the battery charger. Please note the following when using the charger:

Only use a battery charger supplied by Ottobock which has been tested and approved for the respective batteries by Ottobock (observe information on the charger). Failure to do so can result in a battery explosion and possible impairment of health due to contact with battery acid.

- The information on the nameplate of the charger must match the country-specific voltage of the respective power supply network.
- Only use the battery charger within the specified ranges of temperature and humidity.
- Place the rubber feet of the battery charger on a level surface.
- When setting the charger up close to a window, protect it against direct sunlight.
- Keep the charger from overheating. The vent openings on the back of the charger must not be covered.
- Turn off the control unit during the charging process so that the entire charging current flows to the battery.
- Avoid dust and dirt. Only clean with a dry cloth.

The charger is designed for the 12 V gel batteries installed in the C1000 DS power wheelchair. A characteristic curve designed for the installed battery type is programmed into the charger at the factory. If the battery charger is to be used with another power wheelchair or if new batteries are installed, the setting has to be checked.

Please see the instructions for use supplied with the battery charger for further details on use and on the LED indicators.

To charge the batteries, proceed as follows:

- Switch off the control unit of the C1000 DS power wheelchair.
- 2. Open the cover on the charging receptacle. The charging receptacle is integrated into the right-hand steering motor moulding (see Figure 25).
- 3. Connect the plug of the battery charger to the charging receptacle on the power wheelchair.
- Connect the battery charger to a wall socket and turn it on. Charging starts automatically. The battery capacity can be monitored using the LEDs on the charger.

After the charging process is complete, turn off the charger as follows:

- 5. First disconnect the battery charger from the mains supply. Then disconnect the charging plug.
- 6. Close the charging receptacle cover to protect it against moisture.
- 7. Turn on the control unit. The power wheelchair is now ready for operation.

Figure 25 Charging receptacle

- 1 Charging receptacle (Figure shows version for batteries with 63 Ah)
- 2 Battery charger plug

The battery charger indicates the following states:

| Display | | Function |
|---------------------------|-----|----------------------------|
| Green LED flashing | 1 x | Low-maintenance battery |
| | 2 x | Maintenance-free battery |
| Steady | | Battery is charged to 80 % |
| Yellow LED is illuminated | | Battery is charging |
| Green LED is illuminated | | Battery is fully charged |

| Display | Function |
|------------------------|--|
| Red LED is illuminated | Battery is defective or missing or charging time has been exceeded |

Table 5 Battery charging statuses during the charging process

If none of the LEDs are lit, there is no power supply.

6.6.3 Removing the battery cover

Proceed as follows:

- 1. Loosen the two knurled head screws on the left and right of the drive unit sustainer.
- 2. Remove the battery cover by pulling it forward.

Figure 26 Removing the battery cover

6.6.4 Replacing the batteries

⚠ WARNING

Risk of tipping due to changing the battery position. Ensure that the batteries are up against the stop during installation and that there is no play between the batteries and track when installing the track.

NOTICE

Unauthorised battery replacement. Batteries may only be replaced by a specialist dealer.

Proceed as follows:

- 1. Turn off the control unit.
- 2. Remove the fuse (see Section 9.2).
- 3. Remove the battery cover (see Section 6.6.3)
- 4. Unlock the release bolt and remove the track.

Figure 27 Unlocking the battery drawer and removing the track

- 5. Use a ring or open-end wrench size 11 to disconnect all cables from the batteries.
- 6. Remove the batteries by the handles.

Figure 28 Removing the batteries by the handles

- 7. Put the new batteries in place. Ensure that the batteries are up against the stop.
- 8. Use a ring or open-end wrench size 11 to connect all cables to the batteries. Observe correct polarity (black cable = connection to negative pole, red cable = connection to the positive pole).
- Reinstall the track. Ensure that there is no play between the batteries and track when installing the track.

Figure 29 Battery stop

1 Stop, rear position

INFORMATION

When installing the track, ensure the release bolt engages correctly. It has to lock into the corresponding recess on the drive unit sustainer.

Figure 30 Release bolt engaged

6.6.5 Changing the battery installation position

⚠ WARNING

Risk of tipping due to changes in the centre of gravity. Changing the battery position also changes the centre of gravity of the C1000 DS. There is a risk that the steering casters may loose ground contact and the front anti-tipper may engage.

Test the effects of the centre of gravity shift on the behaviour of the power wheelchair with the support of an attendant.

⚠ WARNING

Risk of tipping due to changing the battery position. Ensure that the batteries are up against the stop during installation and that there is no play between the batteries and track when installing the track.

NOTICE

Unauthorised battery replacement. Changing the battery installation location may only be completed by a specialist dealer.

By default, the batteries in the C1000 DS power wheelchair are installed in the far forward position. Where required according to the indications, e.g. dwarfism, the installation position of the batteries in the battery drawer must be changed.

Changing the installation position of the batteries shifts the centre of gravity of the power wheelchair. The driving characteristics of the C1000 DS power wheelchair were tested with all three possible battery positions and conform to the test standard requirements in all cases.

Proceed as follows:

- Turn off the control unit.
- 2. Remove the fuse (see Section 9.2).

- 3. Remove the batteries (see Section 6.6.4)
- 4. Loosen the two Allen head screws on the stop using an Allen wrench size 4 and remove the stop.
- Move the stop according to the mounting bores on the bottom of the battery drawer and reattach the stop using the two Allen head screws.
- 6. Put the batteries back into place.
- 7. Reinstall the batteries (see Section 6.6.4).

Figure 31 Battery stop

- I Mounting holes
- 2 Allen head screws
- 3 Stop

7 Accessories

⚠ CAUTION

Risk of accidents due to unsecured screw connections. When screw connections with thread locks are loosened, they must be replaced with new ones or secured with a medium-strength liquid thread lock substance (e.g. Loctite 241®. After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened. Observe the applicable torque specifications.

INFORMATION

Only original options provided by the manufacturer may be used. The optional components may only be mounted as described here. Failure to comply will void the warranty.

INFORMATION

All available optional add-on components are contained on the order form and in the wheelchair accessories catalogue.

The C1000 DS power wheelchair has been designed as a modular system. Certain component groups can be exchanged and other accessories can be added.

The enAble50 control unit facilitates the use of numerous additional electrical functions as well as customised input devices.

The complete range of options is listed on the order form and in the wheelchair accessories catalogue.

7.1 Seat options

▲ WARNING

Risk of pinching or crushing in the electric seat adjustment range. When activating the fixture to adjust the seat height, back angle and seat tilt, there are pinch and shear points in the area between the seat frame and wheelchair frame due to the construction. All attendants must be informed about this. In order to prevent injuries, ensure that body parts such as the hand and feet are always kept out of the danger zone, that no interfering objects such as clothing or obstacles are in the danger zone and that no unauthorised persons are present there.

⚠ CAUTION

Risk of injury if the power wheelchair tips over. The seat tilt must be in the home position, the fixture to adjust the seat height must be in the lowest position and the back rest must be upright before negotiating uphill slopes and driving up low curbs.

When negotiating downhill slopes or driving down low curbs, moving the seat tilt back slightly and approaching at low speed (max. 3 km/h) is recommended.

INFORMATION

The seat function actuators are not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time). General recommendation: At maximum load capacity, 10 seconds of operation should be followed by 90 seconds of idle time. The electric seat functions are considered independently of the driving function for this purpose.

INFORMATION

When getting into or out of the power wheelchair, bring the seat to a level position.

INFORMATION

For safe operation, please also observe the information in Section 2.3 "General Safety Information" and Section 2.5 "Safety Requirements for Operation".

7.1.1 Seat height adjustment

▲ WARNING

Risk of injury if the power wheelchair tips over. Creep speed is automatically activated when the seat height adjustment function is operated. If this is not case, it represents a malfunction and the seat height adjustment function must not be used. Contact a specialist dealer as soon as possible to eliminate the problem.

⚠ CAUTION

Risk of injury if the power wheelchair tips over. When the seat height has been raised, the power wheelchair may tip over and the user may fall out. The following safety instructions must be followed under all circumstances:

■ When the seat height adjustment fixture is activated, only cover short distance and / or manoeuvre.

- The maximum load capacity of the seat with the seat height adjustment function is 130 kg (optional 200 kg).
- The seat height adjustment function may only be used with the back rest in the upright position.
- Apply the lap belt if possible.
- Use the seat height adjustment function only on level ground.
- Do not modify the installation position of the seat height adjustment fixture.

⚠ CAUTION

Risk of falling. Do not lean the upper part of the body beyond the seat area while the seat is raised.

⚠ CAUTION

Risk of accidents due to a restricted field of vision. Observe the restricted field of vision when the seat is raised and when operating the seat with integrated height adjustment.

NOTICE

Risk of damage due to overloading/excessive strain. The maximum load capacity for the standard version of the C1000 DS power wheelchair is 140 kg. It is possible to

have the power wheelchair retrofitted for a maximum load capacity of 200 kg. If the power wheelchair is equipped with a seat height adjustment unit, the maximum load capacity is 130 kg, optional 200 kg.

NOTICE

Risk of damage due to improper maintenance. Inspect the seat height adjustment fixture for visible signs of damage at least once a month and ensure all screw connections are tight. Also maintain correct air pressure in the tyres.

NOTICE

Risk of transportation damage. Always lower the seat height adjustment fixture fully for loading or transportation.

INFORMATION

A sufficient supply of power must be maintained in order to guarantee proper operation of the seat height adjustment function.

INFORMATION

The seat height adjustment fixture is not intended for continuous operation but only for a 10 % load. For an operating time of 1 hour, the seat height adjustment function can be operated for a maximum of 6 minutes.

The C1000 DS power wheelchair can be equipped with a seat height adjustment fixture as an option. The seat height adjustment function raises the seat bottom above the drive motor by up to 40 cm. It is controlled by the joystick in "additional function" mode or by the push button module.

The driving function remains available even when the seat is raised. As soon as the seat leaves the lowest position, the speed is reduced (creep speed) because the power wheelchair is less tip-resistant. This is indicated by a snail symbol on the control panel or on the LCD monitor if a separate LCD monitor is used.

Figure 32 Creep speed indicator on the control unit display / LCD display

Figure 33 Seat with activated seat height adjustment function

7.1.2 Electric seat tilt

INFORMATION

Please observe the information at the beginning of Section 7.1 "Seat options" for safe operation.

The C1000 DS may be equipped with an optional seat tilt function. The electric seat tilt adjustment lets you bring the seat into a tilted position of up to 30°, e.g. for pressure relief. It is controlled by the joystick in "additional function" mode or by the push button module. The seat can be tilted backwards continuously. The power wheelchair has an integrated feature that shifts the centre of gravity and serves to improve stability against tipping.

Figure 34 Electric seat tilt

7.1.3 Electric back angle adjustment

INFORMATION

Please observe the information at the beginning of Section 7.1 "Seat options" for safe operation.

The back rest can be equipped with an electric back angle adjustment option (see Figure 35). It is controlled by the joystick in "additional function" mode or by the push button module.

In the sitting position, the back rest can be inclined backwards continuously by up to 30°.

Figure 35 Electric back angle adjustment

7.1.4 Recaro® seat

INFORMATION

The conversion of the C1000 DS power wheelchair to a Recaro® seat may only be completed by authorised personnel.

The C1000 DS power wheelchair can be equipped with various optional Recaro® seat models.

Figure 36 C1000 DS power wheelchair with Recaro® seat

To adjust the back angle of a Recaro® seat model, use the rotary handle on the right or left side of the back rest.

A release handle is used to fold down the back rest. It is located on the side of the back rest. By pulling the release handle upward, the back rest can be folded down to the front. The back rest automatically locks in place when it is brought into its initial position. Verify that the back rest is securely locked in place every time it has been unlocked.

The Recaro® seat can be equipped with an electric back angle adjustment option. This function is activated using the control unit in the "additional function" mode.

With the Recaro® LT back rest, the lateral guide in the lumbar section can be individually adjusted by turning the hand wheel. Turning the wheel to the front evenly moves the two lateral guides closer to each other. Turning the wheel to the rear adjusts the two lateral guides away from each other.

Seat bottom

The type X and W seat bottoms have an extendable pad at the front of the seat bottom. The lever used to adjust this pad is located below the seat bottom. Pull the lever up to shift the pad into the desired position. The pad locks into place as soon as the adjustment lever is released. The adjustment lever must audibly lock in place and return to its initial position.

Headrest

The headrest is adjustable in both height and angle.

A second person is needed to remove the headrest. The push buttons must be detected under the cover by feel and pressed simultaneously. Now the second person can pull the headrest up and out.

Figure 37 Illustration of the push buttons

Removing the Recaro® seat from the chassis

- 1. Pull the release strap on the right and left under the seat bottom forward. This loosens the right and left bolts that serve to lock the seat to the frame.
- To release the posterior seat adapter from the retainer bushings, slightly tilt the seat to the rear and slide back a little.
- Remove the Recaro[®] seat.

Attaching the Recaro® seat to the chassis

⚠ CAUTION

Risk of injury if the locking bolts are not engaged. When mounting the Recaro® seat, the seat adapter must be securely attached to the retainer bushings. The anterior locking bolts have to snap into place up to the key ring.

- 1. Place the posterior edge of the Recaro® seat on the end of the seat frame.
- 2. Slide the seat forward until the posterior seat adapter hangs in the retainer bushings.
- 3. Tilt the seat to the front until the locking bolts engage with the seat frame up to the key ring.

7.1.5 Contour seat

INFORMATION

The conversion of the C1000 DS power wheelchair to a contour seat may only be completed by authorised personnel.

The C1000 DS power wheelchair can be equipped with an optional contour seat. Various sizes are available. This results in enhanced sitting comfort and additional adaptation possibilities.

Figure 38 Contour seat

7.2 Push button module

⚠ CAUTION

Risk of accidents due to uncontrolled driving characteristics. The push button module is an option for the use of the enAble50 system. Improper adaptation of the push button module to a power wheelchair with components outside the Ottobock modular system or subsequent modifications by the user are not permitted.

The push button module lets the user select additional electrical functions directly during normal operation.

Accessible functions

Depending on the design of the power wheelchair and programming of the push button module, up to 5 electrical functions can be controlled during normal operation (Figure 39, item 1). Examples include:

- Seat tilt
- Electric back angle adjustment
- Lift seat function
- Right / left footrest; combination of both footrests
- Special functions (combinations) => symbols S1 S5

The symbols used on the push button module correspond to the symbols used on the LCD monitor (see 7.9.1 "Separate LCD monitor with infrared).

Push button functions

The "M" button (= mode, Figure 39, item 2) is used to advance the individual functions (1-2-3-4-5-1-2-...).

The blue LED (Figure 39, item 3) indicates the selected function.

The selected function is executed with the arrow keys (Figure 39, item 4) (function up / down).

INFORMATION

The function of the arrow keys can be adapted according to the wishes of the user (up arrow = function up or function down – depending on programming). Programming may only be completed by authorised personnel.

Three freely positioned Buddy buttons may also be connected to the push button module as options (Figure 39, item 5). They assume the function of the mode and up / down arrow keys (Figure 39, item 2/4). Symbols on the back of the push button module indicate the corresponding functions.


Figure 39 Push button module

- 1 Available electrical functions
- 2 "Function selection" mode button
- 3 "Selected function" LED indicator
- 4 Up / Down arrow keys (see information box)
- 5 Stereo plug for the connection of Buddy buttons

7.3 2x adapter cable for Piko or Buddy button

The C1000 DS power wheelchair may be equipped with a 2x adapter cable. The respective Piko or Buddy button that is connected may be used to control electrical seat functions. The functions to be controlled are selected when the C1000 DS power wheelchair is ordered. Your specialist dealer can change the assignment of the button using a handheld programming device or programming station.

The buttons are equipped with a so-called toggle function. This means that the selected function is activated as long as the button is pressed.

Pressing the button again activates the selected function in the opposite direction. For example, when the function for the electric back rest is activated, the back rest moves back as long as the button is pushed the first time the button is activated. Activating this function again causes the back rest to move forward as long as the button is pushed.

7.4 Mechanically elevating footrests

The C1000 DS power wheelchair may be equipped with mechanically elevating footrests.

To elevate the footrest, proceed as follows:

- Activate the release lever on the footrest (the gas compression spring is activated, see arrow in Figure 40).
- 2. Move the footrest into the desired position.
- 3. Allow the release lever to engage (the gas compression spring is deactivated).

Figure 40 Mechanically elevated footrest

7.5 Electrically adjustable footrests

INFORMATION

The seat function actuators are not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time). General recommendation: At maximum load capacity, 10 seconds of operation should be followed by 90 seconds of idle time. The electric seat functions are considered independently of the driving function for this purpose.

INFORMATION

The conversion of the C1000 DS power wheelchair to electrically adjustable footrests may only be completed by authorised personnel.

To avoid permanent pressure loads or to ensure anti-shock positioning, the C1000 DS power wheelchair offers the option of mounting electrically adjustable footrests. With the control unit in "additional function" mode, the right and left footrest can be operated separately or simultaneously. To adjust them, use the joystick.

The amount of space available for getting into or out of the wheelchair can be increased by flipping the footplates up.

Figure 41 Elevated footrest

7.6 Control panel holder

7.6.1 Swing-away control panel holder with removable control panel

The removable control panel can be pulled up and off the control panel holder.

To allow the user to drive the power wheelchair closer to an object or under the edge of a table, the control panel can be swung to the side with a special control panel holder (see fig. 42).

- 1. Apply slight pressure to push the control panel holder to the side and release the pivot element.
- 2. Swing the control panel holder to the side.
- 3. When bringing the control panel holder back to its original position, the pivot element automatically locks into place.

Figure 42 Control panel holder, swing-away

7.6.2 Height-adjustable control panel holder

An optional special control panel holder can be mounted on the power wheelchair, allowing the height of the control panel to be adjusted (see Figure 43).

- 1. Loosen the two mounting screws.
- Adjust the height of the control panel by moving the control panel up or down on the track.
- 3. Retighten the two mounting screws.

Figure 43 Control panel holder, height-adjustable

7.7 Special controls

Various special controls such as the sip and puff control, chin control or mini joystick can be added to the C1000 DS power wheelchair.

INFORMATION

Detailed information can be found in the instructions for use "Special Controls", reference number 647G646*.

7.8 Attendant control

A separate control panel can be used on the wheelchair for operation by the attendant. It is installed on the push handle and is height adjustable / removable.

Fig. 44 Attendant control

7.8.1 Functional overview

The attendant control can be used to operate the driving function and the electric seat functions.

The module can be connected to the control unit in combination with the hand control device or separately.


Fig. 45 Attendant control overview

- 1 On / Off button
- 2 Battery LED
- 3 Mode button
- 4 Mode LFD

- 5 Display of the selected electric seat function
- 6 Joystick

On/Off button (item 1)

The button is used to turn the power wheelchair on, to activate the drive-away lock (also see chapter 6.3.7), and to turn the wheelchair off.

The control panel displays the following after the attendant control is turned on:

| Display symbol | Information |
|----------------|-----------------------------|
| | Attendant control activated |

Tab. 5 "Attendant control activated" display

Battery LED (item 2)

The battery LED lights up when the unit is on. Colour and flash codes provide information on the current battery capacity status:

| LED | Status | |
|-------------------------|----------------------------|--|
| Green is lit | Battery capacity > 70% | |
| Orange is lit | Battery capacity 30%-70% | |
| Red is lit | Battery capacity < 30% | |
| Pad is flashing | Battery deep discharge - | |
| Red is flashing | charge immediately | |
| | Battery over-voltage (e.g. | |
| Green is flashing | after driving downhill) - | |
| | Continue driving slowly | |
| Red/Orange/Green flash- | Charging process / drive | |
| ing alternately | away lock | |

Tab. 6 "Battery Capacity" indicator

Mode button (item 3)

Switch between speed level 1 and 2 by pressing the button briefly.

Push and hold the button (at least 2 seconds) to switch the control unit to the electric seat functions.

Mode LED (item 4)

The Mode LED indicates the speed level status and system faults using 3 colour LEDs and various flash codes:

| LED | Status |
|----------------|----------------------------|
| Green is lit | Speed level 1 |
| Orange is lit | Speed level 2 |
| Red – flashing | Fault (see chapter 8.2.2) |
| 1x/2x/3x/4x/5x | r duit (see chapter 0.2.2) |
| Off | LED display switched to |
| <u> </u> | seat functions |

Tab. 7 Display of the speed levels / system faults

Display of the selected electric seat function (item 5)

The electric seat functions are accessed by pressing and holding the "M" button (approx. 2 seconds).

The currently selected seat function is indicated on the control panel display (see chapter 6.3.6) and by the following LEDs on the attendant control:

| LED indicator | Information | |
|---------------|---|--|
| | Electric back angle adjust- ment | |
| | Electric seat tilt | |
| | Electric seat height adjust- ment | |
| | Coupled electric seat adjustment (backrest and seat tilt) | |

Tab. 8 Display of electric seat functions

Moving the joystick right switches between the different seat functions. The selected function can be extended or retracted by moving the joystick forward or backward (also see chapter 6.3.6).

The corresponding LED flashes if an error occurs. See chapter 8.2.2 for the meaning of the fault notifications.

Joystick (item 6)

The joystick is used for driving. The further it is moved away from the centre position, the faster the power wheelchair will drive in this direction.

If the seat functions are accessed with the Mode button, then the joystick is used to select and activate the seat functions (see item 5 and chapter 6.3.6).

7.9 Additional indicator and control elements

7.9.1 Separate LCD monitor with infrared

⚠ CAUTION

Risk of injury due to faulty installation of third-party options. Ottobock includes an interface for alternate input devices (environment and special controls) with the LCD monitor. Ottobock does not assume any liability for combinations involving third-party special controls.

The LCD monitor with infrared allows all types of infrared devices to be controlled using the joystick.

Standard household and HiFi devices can be saved in the control unit using a learning mode. Switches or dimmers are offered as accessories. The mouse emulator is another

available infrared module. It allows joystick signals to be used in order to control a computer mouse.

The input devices of the special controls (e. g. sip and puff control, chin control) can also be used for infrared functions.

INFORMATION

Detailed information is contained in the "Special Controls" instructions for use, reference number 647G636*.

Figure 46 LCD module with infrared

| Display | Function |
|--------------|---|
| 14,2 km/h | Driving menu with speed level and battery capacity, left direction indicator light activated |
| WARNING EXX | Driving menu with speed level and low battery capacity |
| ERROR 111 | Charging process with drive stop |
| 3,5 g | Creep speed |

| Display | Function |
|-------------|--|
| | Electric seat height adjust- ment (optional) |
| | Electric back adjustment (optional) |
| | Electric seat tilt (optional) |
| (B) (S) (S) | Special functions, e.g. coupled electric seat adjustment (back adjustment and seat tilt) |

| Display | Function |
|---------|--|
| | Electric footrest left (optional) |
| | Electric footrest right (optional) |
| | Electric footrests, coupled (optional) |
| | Drive-away lock |

Table 9 Information display on the LCD monitor

7.9.2 External kilometre indicator

The power wheelchair can be equipped with an external odometer.

It is connected to the protective bow for the control panel.

Additional information on using it can be found in the manufacturer's instructions for use.

Figure 47 Kilometre indicator

7.10 Additional options

7.10.1Lap belt

⚠ WARNING

Risk of injury due to improper use of the lap belt. The lap belt helps to additionally stabilise the person sitting in the power wheelchair. The lap belt must never be used as part of a restraint system for transportation of the person in a motor vehicle for the disabled.

INFORMATION

For the personal safety of the user, wearing the lap belt at all times is recommended. The lap belt helps to additionally stabilise the person sitting in the power wheelchair.

INFORMATION

If the C1000 DS power wheelchair is designed for speeds of 10 km/h and up, the lap belt is included as standard equipment.

The C1000 DS power wheelchair can be equipped with an optional lap belt. The lap belt serves exclusively as a safety belt when driving the C1000 DS power wheelchair.

Applying the lap belt

- 1. Push the two halves of the buckle together so they engage. The belt buckle must snap in place audibly.
- Verify that the belt has locked by trying to pulling it apart. The lap belt should not be too tight on the body. Objects caught under the belt can cause painful pressure sores.

Opening the lap belt

1. Press the red release button. The belt opens.

Changing the belt length

The belt length can be adjusted on both sides.

- 1. Position both buckle halves in the middle of the body.
- After positioning the buckle halves at a right angle in relation to the belt, their position can be varied. Any excessive belt length will be held in place by the plastic slides.

Figure 48 Adjusting the lap belt

Belt buckle

7.10.2Four-way chest strap/chest belt (for Recaro® seats only)

⚠ WARNING

Risk of injury due to improper use of the four-way chest strap. The four-way chest strap / chest belt helps to additionally stabilise the person sitting in the C1000 DS power wheelchair. The four-way chest strap / chest belt must never be used as part of a restraint system for transportation of the person in a motor vehicle for the disabled.

When the C1000 DS power wheelchair is equipped with a Recaro® seat, it may be retrofitted with a four-way chest strap or chest belt.

The four-way chest strap makes it possible to restrain the patient in the C1000 DS power wheelchair. The lower strap should be placed between the pelvic bone and thigh, and should not be too tight. The belt buckle is positioned at the front of the body, in the centre. The shoulder straps should be positioned evenly on both sides on the thorax, and should not be too tight.

Applying the four-way chest strap / chest belt

- To close the four-way chest strap, connect all parts of the belt buckle. The belt buckle must snap in place audibly.
- Verify that the belt has locked by trying to pulling it apart. The four-way chest strap / chest belt should not be too tight on the body. Objects caught under the belt can cause painful pressure sores.

Changing the belt length

To increase the length of the belt, turn the buckle by 90° and pull. To reduce the length of the belt, pull the overhanging belt end.

7.10.3Headrest mounting adapter (standard seat, contour seat)

The C1000 DS power wheelchair can be equipped with an optional headrest. An adapter is available for the headrest installation kit.

Figure 49 Adapter for headrest installation kit

7.10.4Other optional add-on components

- Puncture-proof tyres: Solid rubber tyres
- Seating shell interface: For the adaptation of special seating shells. Control panel holders for seating shells are also available.

- Arm rest accessories: Special adapter for the arm rests in our wheelchair accessories catalogue
- Joystick accessories: Fork for tetraplegics, STICK S80, soft ball, golf ball and flexible control stick shaft
- Protective bow for control panel: Metal bow to protect the control panel against impacts
- Crutch holder
- Folding rear view mirror
- Tray
- Attachable table tray

INFORMATION

These and other optional add-on components are included on the order form and in the wheelchair accessories catalogue.

8 Malfunctions / troubleshooting

⚠ CAUTION

Risk of accidents and injury due to uncontrolled movements of the power wheelchair. Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunctions. In this case, please contact your authorised dealer immediately.

INFORMATION

During the use of the power wheelchair, electrical discharges may occur which are caused by factors such as friction (high voltage with low current; discharge via the user). However, these do not represent a health hazard.

This reduction in comfort – under consideration of the environmental conditions the user is exposed to – can be eliminated by custom fabrication measures (installing a mechanical discharge contact / ground strap to the wheelchair frame).

Electrostatic discharge may also occur if the power wheelchair is equipped with puncture-proof tyres. Retrofitting the wheelchair with pneumatic tyres can correct this problem.

INFORMATION

Should you encounter problems while troubleshooting or if you do not manage to completely eliminate a problem by following the measures described here, please contact your specialist dealer.

Malfunctions are indicated on the control panel LCD monitor. Table 6 describes the individual error codes with the corresponding problem sources as well as the possible causes and resolutions.

If malfunctions cannot be fully corrected using the procedures described here, an authorized dealer has the ability to read the exact error code with a handheld programming device and to carry out a targeted system analysis.

All problems that have ever occurred are saved in a list and can be retrieved, e.g. in case of a general overhaul of the power wheelchair. The saved data can be used to determine future service and maintenance intervals, for example.

8.1 Warning

A warning indicates a status or malfunction of one or several components of the power wheelchair. The function of components that have no errors is not restricted. For example, if the connection between the controller and seat motor malfunctions, this error will only be indicated if the user attempts to activate the motor. However, the driving function is still available.

8.2 Error

⚠ CAUTION

Risk of injury if the power wheelchair stops abruptly. In the event of communication problems in the bus system of the control unit, the system triggers an emergency stop to prevent any uncontrolled functions. Depending on the source of the malfunction, it is possible to drive the wheel-chair out of a danger zone such as road traffic after switching the control unit on again.

If the driving function is still not available after switching the control unit on again, the brake must be deactivated (see Section 6.5) to switch over to pushing mode. After the power wheelchair has triggered an emergency stop, contact a specialist dealer as soon as possible!

⚠ CAUTION

Risk of accidents due to uncontrolled driving characteristics. Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunctions. In this case, please contact your authorised specialist dealer immediately.

An error affects one or several functions of the power wheelchair. The system is not fully operational until the error has been corrected.

8.2.1 Error Indication: Control Panel / LCD Monitor

| Display symbol(s) | LCD monitor display | Warning / error | Cause | Possible corrective action |
|-------------------|---------------------|-----------------------------------|---|---|
| | WARNING XXX | Controller temperature warning | Overheating due to | Cool down phase |
| | WARNING XXX | Motor temperature warning | excessive load | Cool down phase |
| | ERROR MX | Joystick warning | Joystick not in zero position when switching on | Bring joystick to zero position before switching on |

| Display symbol(s) | LCD monitor display | Warning / error | Cause | Possible corrective action |
|-------------------|-----------------------|---|---|--|
| × | ERROR XII | Hand control device fault | Defective joystick | Contact enceiglist degler |
| | ERROR 1920 1420 | Controller error | Defective controller | Contact specialist dealer |
| | ERROR MAX | Communication error (alternating flashing signal) | Faulty connection between the hand control device and controller. Defective cabling, software or hardware | Check cabling / plug connections Contact specialist dealer |
| | WARNING LEXT | Battery under voltage | Battery deep discharge | Charge as soon as possible |
| | WARNING AXX 14.25 | Battery over voltage | Voltage too high (after full charge and driving downhill) | Continue driving slowly |

| Display symbol(s) | LCD monitor display | Warning / error | Cause | Possible corrective action |
|---|---------------------|--------------------------------------|---|--|
| | ERROR XXX | Back angle adjustment motor fault | Faulty cabling or plug connection; | Check cabling / plug connections |
| | ERROR XXX | Seat tilt motor fault | defective actuator | Contact specialist dealer |
| | ERROR XXX | Seat height adjustment motor fault | | |
| المالية | ERROR XXX | Electric footrest motor fault | Faulty cabling or plug connection; defective actuator | Check cabling / plug connections Contact specialist dealer |
| | ERROR XXX | Drive motor fault | | |

| Display symbol(s) | LCD monitor display | Warning / error | Cause | Possible corrective action |
|-------------------|----------------------|------------------|---|---|
| | ERROR XII | Wheel lock fault | Wheel lock release open Defective wheel lock | Close wheel lock release Check wheel lock (e.g. Bowden cable) |
| | STOP ERROR EXX | Emergency stop | Severe fault caused by controller, hand control device and / or drive motor malfunction | Contact specialist dealer |

Table 10 Status and error messages

8.2.2 Error Indication: Attendant Control

| LED indicator | Warning/Error | Cause | Possible corrective action | |
|--------------------------------|---------------------------|---|----------------------------|--|
| Red is flashing | Battery under voltage | Battery deep discharge | Charge as soon as possible | |
| Green is flashing | Battery over voltage | Voltage too high (after full charge and driving downhill) | Continue driving slowly | |
| Mode LED (see Fig. 45, item 4) | | | | |
| Red is flashing 1 x | Controller fault | Defective controller | Contact specialist dealer | |
| Red is flashing 2 x | Hand control device fault | Defective joystick | Contact specialist dealer | |

| Red is flashing 3 x | Communication error | Faulty connection between the hand control device and controller; | Check cabling / plug connections; |
|--|--|---|---|
| | | Defective cabling, software or hardware | Contact specialist dealer |
| Red is flashing 4 x | Joystick warning | Joystick is not in zero position when the unit is turned on | Move the joystick to the zero position before turning the unit on |
| | Brake fault | Open brake release; | Close brake release; |
| Pod is flashing 5 v | | Defective brake | Check brake (e.g. Bowden cable) |
| Red is flashing 5 x | Drive motor fault | Faulty cabling or plug connection; | Check cabling / plug connections; |
| | | defective actuator | Contact specialist dealer |
| Seat function LEDs (see Fig. 45, item 5) | | | |
| Backrest LED is flashing | Back angle adjustment motor fault | | |
| Seat bottom LED is flashing | Seat tilt motor fault | Faulty cabling or plug contact; | Check cabling / plug connections; |
| 2 LEDs below seat bottom are flashing | Seat height adjustment motor fault | Defective actuator | Contact specialist dealer |
| Backrest / seat bottom LED is flashing | Electric seat function temperature warning | Overheating due to excessive load | Cool down phase |

Table 11 Attendant control error messages

8.3 Defect/failure

A defect represents a serious failure of a system component. A failure is the most serious type of malfunction, which causes an immediate emergency stop of the system.

Defects / failures are shown by a continuous sequential light indicator with an audible signal.

After the malfunction has been corrected, the system is reactivated by restarting it.

9 Maintenance, cleaning and disinfection

INFORMATION

A spare parts catalogue for ordering spare parts is available from Ottobock upon request. Only original Ottobock replacement and wear parts must be used. Failure to comply with this requirement will render the warranty null and void.

INFORMATION

Should you encounter problems during maintenance, contact your specialist dealer. The operational safety power wheelchair is to be checked once a year by a specialist dealer.

9.1 Maintenance intervals

The correct functionality of the power wheelchair must be checked prior to each use. The items listed in Table 7 must be checked by the user at the indicated intervals.

| Component | Activity | Daily | Weekly | Monthly |
|--------------|---|-------|--------------------|---------|
| Arm rest and | Fastening screws tightened | | | Х |
| side panel | Arm rest and control panel secured | | Prior to every use | |
| | Check arm rest for damage | | X | |
| Drive wheels | Wheels must rotate freely and without axial run out | | | Х |
| | Check if attachment screws are seated securely | | | X |
| | Check if wheel mounts are seated securely | | | X |
| | Directional stability of the entire power wheelchair | | X | |
| Tyres | Air pressure (printed on the sidewall of the tyre) | | | Х |
| | Sufficient tread depth, at least 1 mm | | | X |
| | Check for damage | | | X |
| Batteries | Maintenance-free | | | |
| Lighting | Visually check for damage on the housing | | Х | |
| | Verify electric function | X | | |
| Electronics | Control unit / push button module free of error messages | | Prior to every use | |
| | Charger does not indicate any error messages on the LCD monitor | | X | |
| | Check plug connections | | | Х |

| Component | Activity | Daily | Weekly | Monthly |
|------------------|---|-------|--------|---------|
| Wheel lock | Activate wheel lock lever while control unit is switched on | X | | |
| | Brake function active with wheel lock engaged | | | X |
| Footrest | Check ratchet mechanism for functionality | | | |
| | and secure seating | | | X |
| | Check footrests for damage | | | X |
| Electric seat | Visual inspection of all moving components, especially | | | |
| adjustments | cabling - check for damage | | | X |
| | Check whether screw connections are tight | | | X |
| Steering casters | Verify function of the steering motors | | | Х |
| | Unrestricted steering angle | | | X |
| Padding and | Proper condition of padding | | | Х |
| belts | No wear on the seat belts | | | X |
| | Check belt buckle for functionality | | X | |
| Seat attachment | Check if attachment screws are seated securely | | | Х |
| | Check seat lock of Recaro® seat | | | X |

Table 12 Maintenance measures and intervals

INFORMATION

When the joystick is activated while the brake is deactivated, the control unit emits an error signal on the control panel. If this is not the case, there is a malfunction that must be corrected immediately by a specialist dealer.

9.2 Changing the fuse

The 80 A fuse is located in a fuse holder at the rear end of the battery case (see Figure 49).

- 1. Open the cap of the fuse holder.
- Pull out the fuse.
- 3. Insert a new fuse into the holder. Make sure that the fuse is pressed into the centre of the spring contacts and that it is not positioned at an angle.
- 4. Close the cap until you feel it snap back into place.


Figure 49 Fuse holder

- 1 Cover open
- 2 Fuse, inserted

9.3 Wheel replacement

NOTICE

Risk of damage due to uncontrolled movements. When raising the power wheelchair, e.g. with a jack, secure it against slipping and tipping sideways by using a suitable block under the drive unit sustainer.

NOTICE

Damage to tyres due to excessive tyre pressure. Do not exceed the tyre pressure specified in Section 10 "Technical Data".

NOTICE

Damage to the drive unit sustainer. When raising the power wheelchair, e.g. with a jack, avoid mechanical damage to the drive unit sustainer such as surface scratches. Apply the jack to the bottom of the drive unit sustainer.

INFORMATION

Direct sunlight (UV light) causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread.

INFORMATION

Avoid parking the wheelchair outdoors whenever possible. Regardless of wear and tear, replace the tyres every 2 years. If the wheelchair is parked for an extended period of time or the tyres overheat (e.g. in the vicinity of radiators or in case of exposure to strong sunlight behind glass), the tyres may become permanently deformed. Therefore, sufficient

clearance should always be maintained from sources of heat, and the wheelchair should be moved frequently or placed on blocks for storage.

The C1000 DS power wheelchair is propelled by two 14" drive wheels and steered by means of two 10" steering casters.

Figure 50 Jack

9.3.1 Steering caster replacement

Proceed as follows:

- 1. Turn off the control unit.
- 2. Raise the power wheelchair, e.g. using a jack, so that the wheel being replaced can turn freely.

- 3. Remove the cover cap with a flathead screwdriver.
- 4. Loosen the hexagon head screw with a socket wrench size 17.
- 5. Pull the steering caster forward off the wheel hub.
- 6. Push the new steering caster onto the wheel hub.
- 7. Retighten the hexagon head screw.
- 8. Check proper tightness with a torque wrench.
- 9. Reapply the cover cap.

INFORMATION

When reassembling the wheel, tighten the screw to a torque of **25 Nm**.

Figure 51 Removing the steering caster

9.3.2 Drive wheel replacement

Proceed as follows:

- 1. Turn off the control unit.
- 2. Raise the power wheelchair, e.g. using a jack, so that the wheel being replaced can turn freely.
- Loosen the four Allen head screws using an Allen wrench size 6 and remove them including the spring washers.
- 4. Pull the drive wheel forward off the wheel hub.
- 5. Push the new drive wheel onto the wheel hub.

- 6. Put the spring washers on.
- 7. Alternately tighten the four Allen head screws crosswise.
- 8. Check proper tightness with a torque wrench.

INFORMATION

When reassembling the wheel, tighten the screws to a torque of 25 Nm.

Figure 52 Removing the drive wheel

9.3.3 Replacing the casing and/or tube for pneumatic tyres

The wheel rims are in two parts and can be separated by removing the Allen head screws. Proceed as follows:

- 1. Remove the wheel in question from the power wheelchair (see Section 9.3.1 and 9.3.2).
- 2. For pneumatic tyres: Let all the air out of the tyre.
- Loosen all five Allen head screws that connect the two parts of the rim to each other using an Allen wrench size 6. Remove the screws.
- 4. Using a tyre lever, remove the casing from the rim.
- 5. For pneumatic tyres: Press the valve to the inside through the rim.
- 6. For pneumatic tyres: Pull out the tube.
- 7. For pneumatic tyres: Repair the inner tube with a standard tyre repair kit or replace it with a new one.
- If applicable, replace the defective casing with a new one and install it on the rim.
- Join the two parts of the rim using the five Allen head screws.

- 10. For pneumatic tyres: Inflate the tyre; observe the allowable tyre pressure.
- 11. Reinstall the wheel on the power wheelchair (see Section 9.3.1 and 9.3.2).

Figure 53 Removing the casing

9.4 Changing defective lights

9.4.1 Front lights

INFORMATION

The front lights may only be installed and removed by authorised personnel.

INFORMATION

Only the complete front light can be replaced on each side.

An LED headlight and LED direction indicator light are attached to the right and left splash guards of the drive wheels, respectively.

9.4.2 Rear lights

INFORMATION

The rear lights can be ordered from the specialist dealer.

INFORMATION

Only the complete rear light can be replaced on each side.

Two LED rear lights with integrated LED direction indicator lights are installed in the steering motor mouldings at the rear of the power wheelchair.

To replace the rear lights, proceed as follows:

1. Loosen and remove the Allen head screw using an Allen wrench size 4.

Figure 54 Attaching the rear light to the steering caster fork

2. Pull the mounting plate with the rear light forward.

Figure 55 Removing the rear light

Separate the electric supply line from the connector plug.

Figure 56 Rear light connector plug

- 1 Plug connection from the rear light to the controller
- 4. Remove the rear light.
- 5. Take the new rear light and connect the electric supply line to the controller using the connector plug.
- 6. Position the rear light in the moulding.
- Insert the Allen head screw and tighten it using an Allen wrench size 4.

9.5 Cleaning and Care

NOTICE

Risk of damage to electronics due to water penetration. In order to avoid malfunctions, the electronic components, motors and batteries must not come into contact with water when cleaning the power wheelchair.

NOTICE

Risk of damage to power wheelchair components. Do not use any aggressive cleaners, solvents or hard brushes to clean the power wheelchair. Never use a water jet or high-pressure cleaning apparatus to clean the power wheelchair.

INFORMATION

Prior to disinfection, clean the seat and back upholstery as well as the control panel and arm rests.

The power wheelchair must be cleaned regularly, depending on the amount of use and the degree of soiling.

The control panel, battery charger, arm rests and trim can be cleaned with a damp cloth and mild cleaning solution.

Use a dry brush to clean the seat and back upholstery.

Use a damp plastic brush to clean the wheels and frame.

9.5.1 Disinfection

Wipe all parts of the wheelchair down with disinfectant.

Important information about disinfection

- Water based disinfectants should be used. Observe the usage information provided by the manufacturer.
- Prior to disinfection, clean the seat and back upholstery, seat cushion, control panel, and armrest.

10 Technical data

| Dimensions and weights | | |
|------------------------|---|--|
| Seat width | 38 – 42 cm or 43 – 48 cm | |
| Seat depth | 38 – 46 cm or 42 – 50 cm | |
| Seat height | 45 — 60 cm | |
| Arm rest height | 24 – 36 cm | |
| Arm rest length | 26 cm | |
| Lower leg length | 25 – 34 cm or 35 – 44 cm or 45 – 54 cm | |
| Back height | 45 or 55 cm | |
| Back angle | -9/1/11/21° or 0/10/20/30° | |
| Overall width | 65 cm | |
| Overall height | 103 cm | |
| Overall length | 120 cm | |
| Turning radius: | 80 cm | |
| Tyre size | | |
| Front: | 14" | |
| Rear: | 10" | |

| Tyre pressure | Front: Printed on the sidewall of the tyre |
|--|---|
| Weight empty | 129 kg (depending on options) |
| Shipping weight* | see weight when empty*, of which: |
| | Side panel: < 1 kg (2.2 lbs) |
| * The specified weight varies according to the | Footrest (standard equipment): approx. 1 kg (2.2 lbs) |
| selected options and model. | Manual elevating footrest: 1.8 kg (4 lbs) |
| Max. load capacity (with standard battery): | |
| Without seat height adjustment | 140 kg (optional 200 kg) |
| With seat height adjust- ment | 130 kg (optional 200 kg) |

| Max. load capacity (with optional battery): | |
|---|---|
| Without seat height adjustment | 140 kg / 309 lbs (optionally 200 kg / 440 lbs) |
| With seat height adjust- ment | 130 kg / 287 lbs (optionally 200 kg / 440 lbs) |
| Corrosion protection | Coated frame |
| Electrical equipment | |
| Operating voltage | 24 V |
| Batteries: | |
| Gel batteries (standard) | 2 x 12 V; 63 Ah (C5) |
| Gel batteries (option) | 2 x 12 V; 94 Ah (C5) |
| Control unit: | |
| Model | enAble50 |
| Operating voltage | 24 V DC |
| Max. output current per | |
| motor | 2 x 130 A |

| Push button module: | |
|---|---|
| Model | For enAble50 control unit |
| Protection class | IPX4 |
| Operating and storage temperature range | -20 °C to +80 °C |
| Lighting: | |
| Front direction indicator | 12 V, 1 W |
| Front light | 12 V, 2.4 W |
| Rear direction indicator | 12 V, 1 W |
| Rear light | 12 V, 1 W |
| Fuse | 80 A |
| Battery charger | For more information, see the included battery charger instructions for use |

| Driving data | | |
|--|--|--|
| Speed | 6 km/h, 10 km/h, 14 km/h (depending on the model) | |
| Climbing ability | 17 % | |
| Maximum obstacle height | 10 cm | |
| Approx. driving range | 35 km (optional 50 km) | |
| Operating temperature | -25 °C to +50 °C | |
| Transportation and storage temperature range | -40 °C to +65 °C | |

Table 13 Technical data

11 Disposal

INFORMATION

If a wheelchair is to be disposed of, all components and materials of the power wheelchair must be recycled or disposed of properly.

If the power wheelchair is no longer in use, it must be disposed of properly in accordance with national regulations.

Please return defective batteries to your specialist dealer when buying new ones.

12 Information on re-use

The C1000 DS power wheelchair is suitable for re-use.

Similar to second-hand machines or vehicles, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

Based on market observations and the current state-of-the-art, the manufacturer has calculated that the C1000 DS power wheelchair can be used for a period of **5 years**, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the wheelchair is stored at the dealer or with paying parties are not included in this period. It should be clearly pointed out, however, that the C1000 DS is a reliable product far beyond this defined period of time, provided that it is cared for and maintained properly.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected.

Subsequently, the product must be examined by an authorised specialist to check the condition, wear and possible defects.

Any worn and damaged components as well as components which do not fit or are unsuitable for the new user must be replaced.

The service manual includes a service plan for each model, detailed information and a list of the required tools.

13 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

14 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock Australia Pty. Ltd.
Otto Bock Asia Pacific Ltd.
Otto Bock Korea HealthCare Inc. 137-070 Seoul·Korea T +82 (0) 2 577-3831 · F +82 (0) 2 577-3828 info@ottobockkorea.com
Otto Bock Manufacturing Königsee GmbH